







# Venturi Masks

**Intended use:**

Venturi masks are intended for administration of oxygen when precise control of inspired oxygen concentration is mandatory. During inspiration the patient draws gases both from oxygen flowing into the mask through the tubing, as well as from room air via ports on the sides of the mask. Oxygen concentration delivered is depending on the size of diluter. The masks are connected directly to a compressed air or oxygen supply. The products are intended for single use only.



**Product range:**

	REF*	Item	Size	Diluters	Description	Tube length
	93-120MMEU 93-120MMEUV (Vanilla scented mask)	Venturi Mask Kit with Single Diluters	Adult	Single Diluters: 24%, 28%, 31%, 35%, 40%, 50%, 60%	Adult Venturi Mask Kit, 7 diluters (24%, 28%, 31%, 35%, 40%, 50%, 60%) Elongated oxygen mask with nose clip, 2.1m (7') No-crush oxygen tube, Fits-All connector	2.1m (7')
	93-227MMEU 93-227MMEUV (Vanilla scented mask)	Venturi Mask Kit with Single Diluters	Pediatric	Single Diluters: 24%, 28%, 31%, 35%, 40%, 50%, 60%	Pediatric Venturi Mask Kit, 7 diluters (24%, 28%, 31%, 35%, 40%, 50%, 60%) Elongated oxygen mask with nose clip, 2.1m (7') No-crush oxygen tube, Fits-All connector	2.1m (7')
	93-270MMEU	Venturi Mask Kit with Single Diluters	Infant	Single Diluters: 24%, 28%, 31%, 35%, 40%, 50%, 60%	Infant Venturi Mask Kit, 7 diluters (24%, 28%, 31%, 35%, 40%, 50%, 60%) Oxygen mask with nose clip, 2.1m (7') No-crush oxygen tube, Fits-All connector	2.1m (7')
	93-128MMEU	Venturi Mask Kit with Single Diluters	Adult	Single Diluter: 28%	Adult Venturi Mask Kit, 1 diluter (28%) Elongated oxygen mask with nose clip, 2.1m (7') No-crush oxygen tube, Fits-All connector	2.1m (7')
	93-120AMM	Venturi Mask Kit with Adjustable Diluter	Adult	Adjustable Diluter: 24%, 26%, 28%, 31%, 35%, 40%, 50%	Adult, Venturi Mask Kit, Adjustable diluter (24%, 26%, 28%, 31%, 35%, 40%, 50%) Elongated oxygen mask with nose clip, 2.1m (7') No-crush oxygen tube, Fits-All connector	2.1m (7')
	93-227AMM	Venturi Mask Kit with Adjustable Diluter	Pediatric	Adjustable Diluter: 24%, 26%, 28%, 31%, 35%, 40%, 50%	Pediatric Venturi Mask Kit, Adjustable diluter (24%, 26%, 28%, 31%, 35%, 40%, 50%) Elongated oxygen mask with nose clip, 2.1m (7') No-crush oxygen tube, Fits-All connector	2.1m (7')

\* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

**Material specifications:**

Contains latex? ..... No  
 Contains phthalates? ..... No  
 Contains bisphenols? ..... No  
**Mask** ..... Polyvinyl chloride  
**Elastic strap** ..... Polyester  
**Nose clip** ..... Aluminum  
**Oxygen nipple connector** ..... Polyvinyl chloride  
**Diluter** ..... Polyethylene (fixed)  
 Acrylonitrile butadiene styrene, Stainless steel (adjustable)  
**Oxygen tube** ..... Polyvinyl chloride  
**Tube connector** ..... Polyvinyl chloride

**Certifications:**

**CE-marking?** ..... YES (CE 0120)  
**Product classification** ..... IIA  
**Quality system certified to** ..... EN ISO 13485:2012

**Sterilization:**

**Sterile?** ..... No  
**Shelf life** ..... 5 years

## Packaging specifications:

REF	Primary package	Secondary package/Polybag	Tertiary package/Shipper			Pallet
	Qty, pcs	Qty, pcs	Qty, pcs	Dimensions, mm (LxWxD)	Gross Wt., kg	Qty, pcs
93-120MMEU	1	10	50	600 x 400 x 260	5.76	600
93-120MMEUV	1	10	50	600 x 400 x 260	6.30	600
93-227MMEU	1	10	50	600 x 400 x 260	5.22	600
93-227MMEUV	1	10	50	600 x 400 x 260	5.94	600
93-270MMEU	1	10	50	600 x 400 x 260	5.02	600
93-128MMEU	1	10	50	600 x 400 x 260	2.30	600
93-120AMM	1	10	50	600 x 400 x 260	2.85	600
93-227AMM	1	10	50	600 x 400 x 260	2.45	600

## EAN codes:

REF	Primary package	Secondary package	Tertiary package
93-120MMEU	(01)06945099301246	(01)26945099301240	(01)46945099301244
93-120MMEUV	(01)06945099301413	(01)26945099301417	(01)46945099301411
93-227MMEU	(01)06945099301253	(01)26945099301257	(01)46945099301251
93-227MMEUV	(01)06945099301420	(01)26945099301424	(01)46945099301428
93-270MMEU	(01)06945099301260	(01)26945099301264	(01)46945099301268
93-128MMEU	(01)06945099301277	(01)26945099301271	(01)46945099301275
93-120AMM	(01)06945099301284	(01)26945099301288	(01)46945099301282
93-227AMM	(01)06945099301291	(01)26945099301295	(01)46945099301299

## Packaging material:

Primary package	Poly Bag - Low density polyethylene / Insert - 80g Paper
Secondary package	Poly Bag - Low density polyethylene / Sticker - Coated Paper
Tertiary package	Sticker - Coated Paper / Carton - Corrugated Board

## Storage and handling:

	<b>FRAGILE</b> - Contents of the transport package are fragile and shall therefore be handled with care
	<b>KEEP AWAY FROM RAIN</b> - Transport package shall be kept away from rain
	<b>THIS WAY UP</b> - Indicates correct upright position of the transport package
	<b>DO NOT USE HAND HOOKS</b> - Hooks are prohibited for handling the transport package



**Ningbo Shengyurui Medical Appliances Co., Ltd.**

No. 138 Binhaisi Rd.  
Hangzhou Bay New Zone, 315336, Ningbo City  
Zhejiang Province, P.R. China

Made in China



**Shanghai Intl. Holding Corp. GmbH**

Eiffestrasse 80, 20537 Hamburg  
Germany

Distributed by



ConvaTec Limited  
First Avenue, Deeside Industrial Park  
Deeside, Flintshire, CH5 2NU, UK

[www.convatec.com](http://www.convatec.com)

**For further information on these products, please contact your local ConvaTec sales representative.**

With reservations for misprints and alterations in product specifications. Printed version not controlled.  
All trade marks and logos are the property of a ConvaTec group company.

© 2017 ConvaTec Inc. AP-017338-MM, ver. FEB 2017