


Pediatric Foley catheters

100% Silicone, 2-way, hard valve

Intended use:

The pediatric Foley catheter is a urinary catheter intended for short as well as long term drainage of the bladder. For routine and post-operative drainage. The size is specially designed to fit pediatric patients. The product is intended for single use only.



List of product range:

Pediatric Foley catheter, 100% clear Silicone, Standard, 2-way, hard valve

Balloon size	CH 06	CH 08	CH 10
1.5 ml/cc	MM43150603		
3.0 ml/cc		MM43150803	MM43151003

Material specifications:

Contains latex? No
Contains DEHP? Yes (Sleeve only)
Catheter tube and balloon Silicone
Hard valve Polypropylene
Tip Silicone, Barium sulphate
Sleeve Polyvinyl chloride

Certification:

CE-marking? Yes (CE 0123)
Product classification Class IIB (EU)
Quality system certified to ISO 13485
Global Medical Device Nomenclature (GMDN) 34917

Sterilisation:

Sterile? Yes
Sterilized by Ethylene oxide
Shelf life 5 years

Packaging:

REF*	Shelf box		Shipping carton				MOQ, pcs.
	Qty., pcs	Dimensions, mm	Qty., pcs	Dimensions, mm	Gross Wt., kg	Volume, m ³	
MM43151003	10	500 x 92 x 52	500	580 x 500 x 527	10.00	0.15	500
MM43150603	10	500 x 92 x 52	500	580 x 500 x 527	10.60	0.15	500
MM43150803	10	500 x 92 x 52	500	580 x 500 x 527	12.40	0.15	500

* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Storage and handling:



With reservations for misprints and alterations in product specifications. Printed version not controlled.

Manufacturer: Unomedical Sdn Bhd
 (15290-P)
Address: Bakar Arang Industrial Estate
 08000 Sungai Petani, Kedah, Malaysia

European Representative: Unomedical a/s
Address: Kongevejen 2
 3460 Birkerød, Denmark

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.
The ConvaTec logo is a trademark of ConvaTec Inc.