








Oxygen Tubing

Intended use:

Oxygen tubes are used to connect the oxygen source with the oxygenation device during administration of oxygen. The product is intended for single use only.



Product range:

	REF*	Item	Description	Tube length
	93-402MM	No-crush Oxygen Tube	No-crush Oxygen Tube, 2.1m (7') Fits-All connector	2.1m (7')
	93-404MM	No-crush Oxygen Tube	No-crush Oxygen Tube, 4.2m (14') Fits-All connector	4.2m (14')
	93-407MM	No-crush Oxygen Tube	No-crush Oxygen Tube, 7.6m (25') Fits-All connector	7.6m (25')
	93-408MM	No-crush Oxygen Tube	No-crush Oxygen Tube, 15.2m (50') Fits-All connector	15.2m (50')
	93-414MM	Oxygen Tube	Oxygen Tube, 10m (33') ID 4.75mm, OD 7.14mm No connectors	10m (33')
	93-430MM	Oxygen Tube	Oxygen Tube, 30m (100') ID 4.75mm, OD 7.14mm No connectors	30m (100')
	93-969MM	Aerosol-to-Oxygen Nipple Adapter	Aerosol-to-Oxygen Nipple Adapter 22mm ID circuit adapter	-

* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Material specifications:

Contains latex? No
Contains phthalates? No
Contains bisphenols? No
Oxygen tube. Polyvinyl chloride
Connector Polyvinyl chloride
Oxygen reducer Propene polymer

Certifications:

CE-marking? Yes (CE 0120)
Product classification IIA
Quality system certified to EN ISO 13485:2012

Sterilization:

Sterile? No
Shelf life 5 years

Packaging specifications:

REF	Primary package	Secondary package/Polybag	Tertiary package/Shipper			Pallet
	Qty, pcs	Qty, pcs	Qty, pcs	Dimensions, mm (LxWxD)	Gross Wt., kg	Qty, pcs
93-402MM	1	10	50	350 x 300 x 200	2.49	1,600
93-404MM	1	10	50	400 x 400 x 300	4.10	900
93-407MM	1	5	25	380 x 400 x 260	3.72	450
93-408MM	1	10	20	400 x 400 x 300	5.52	360
93-414MM	1	1	10	400 x 300 x 240	2.60	320
93-430MM	1	1	1	200 x 200 x 120	0.90	192
93-969MM	1	10	50	300 x 200 x 200	0.50	3,200





EAN codes:

REF	Primary package	Secondary package	Tertiary package
93-402MM	(01)06945099313058	(01)26945099313052	(01)46945099313056
93-404MM	(01)06945099313065	(01)26945099313069	(01)46945099313063
93-407MM	(01)06945099313072	(01)26945099313076	(01)46945099313070
93-408MM	(01)06945099313089	(01)26945099313083	(01)46945099313087
93-414MM	(01)06945099313096	(01)26945099313090	(01)46945099313094
93-430MM	(01)06945099313102	(01)26945099313106	(01)46945099313100
93-969MM	(01)06945099318015	(01)26945099318019	(01)46945099618013

Packaging material:

Primary package	Poly Bag - Low density polyethylene / Insert - 80g Paper
Secondary package	Poly Bag - Low density polyethylene / Sticker - Coated Paper
Tertiary package	Sticker - Coated Paper / Carton - Corrugated Board

Storage and handling:

	FRAGILE - Contents of the transport package are fragile and shall therefore be handled with care
	KEEP AWAY FROM RAIN - Transport package shall be kept away from rain
	THIS WAY UP - Indicates correct upright position of the transport package
	DO NOT USE HAND HOOKS - Hooks are prohibited for handling the transport package



Ningbo Shengyurui Medical Appliances Co., Ltd.

No. 138 Binhaisi Rd.
Hangzhou Bay New Zone, 315336, Ningbo City
Zhejiang Province, P.R. China

Made in China



Shanghai Intl. Holding Corp. GmbH

Eiffestrasse 80, 20537 Hamburg
Germany

Distributed by



ConvaTec Limited
First Avenue, Deeside Industrial Park
Deeside, Flintshire, CH5 2NU, UK

www.convatec.com

For further information on these products, please contact your local ConvaTec sales representative.

With reservations for misprints and alterations in product specifications. Printed version not controlled.

All trade marks and logos are the property of a ConvaTec group company.

© 2017 ConvaTec Inc. AP-017305-MM, ver. FEB 2017