


HandyVac™

Wound drainage system

Intended use:

Low vacuum wound drainage system for post-operative wound drainage. The product is intended for single use only.



List of product range:

HandyVac™	Drain length mm	Perforation length mm	CH 08	CH 10	CH 12	CH 14	CH 16	CH 18
			Inner/outer diameter of drain, mm					
			1.6/2.7	2.0/3.3	2.2/4.0	2.4/4.7	3.0/5.3	3.5/6.0
			REF*	REF*	REF*	REF*	REF*	REF*
HandyVac™ 700 ml bag, with bottom outlet and spiral drain"	650	150	4009918x	4010118x	4010318x	4010518x	4017118x	4010718x
HandyVac™ 700 ml bag with bottom outlet, solid trocar and spiral drain	650	150	4012118x	4012318x	4016018x	4012518x	4018418x	4012718x
HandyVac™ 700 ml bag without bottom outlet, with solid trocar and spiral drain, with glued connector	650	150	4022118x	4022218x	4022718x	4022318x	4022818x	4022418x
HandyVac™ 700 ml bag without bottom outlet, with solid trocar and spiral drain	650	150	4023218x	4023318x	4023418x	4023518x	4023618x	4023718x
Extra drainage bag 700 ml with bottom outlet			40067069					
Extra drainage bag 700 ml with bottom outlet, 90 cm tube			40069069					
Extra drainage bag 700 ml without bottom outlet			4020018x					

* Here noted as x, reference codes contain an additional digit (1-5) representing a language code.
 Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Material specifications:

Contains latex? No
 Contains DEHP? Yes
 HandyVac™ Bellows Polyethylene
 Drainage bag Polyvinyl chloride
 O-Ring Silicone
 Other components Polyvinyl chloride, Polyethylene,
 Polycarbonate
 Trocar Steel

Certification:

CE-marking? Yes (CE 0086)
 Product classification Class I (EU) (without drain)
 Class IIA (EU) (with drain)
 Quality system certified to ISO 13485
 Global Medical Device Nomenclature (GMDN). 35824

Sterilisation:

Sterile? Yes
 Sterilized by Ethylene oxide
 Shelf life 5 years

HandyVac™

Wound drainage system

Packaging:

REF*	Shelf box		Shipping carton				MOQ, pcs.
	Qty., pcs	Dimensions, mm	Qty., pcs	Dimensions, mm	Gross Wt., kg	Volume, m ³	
4009718x	*	*	20	392 x 306 x 375	3.00	0.045	100
4026618x	*	*	20	392 x 306 x 375	3.00	0.045	100
4009918x	*	*	20	392 x 306 x 375	3.20	0.045	100
4010118x	*	*	20	392 x 306 x 375	3.20	0.045	100
4010318x	*	*	20	392 x 306 x 375	3.20	0.045	100
4010518x	*	*	20	392 x 306 x 375	3.20	0.045	100
4010718x	*	*	20	392 x 306 x 375	3.50	0.045	100
4016118x	*	*	20	392 x 306 x 375	4.00	0.045	100
4020018x	50	294 x 125 x 360	150	392 x 306 x 375	13.30	0.029	150
4017118x	*	*	20	392 x 306 x 375	3.50	0.045	60
4012118x	*	*	20	392 x 306 x 375	3.50	0.045	100
4012318x	*	*	20	392 x 306 x 375	3.60	0.045	100
4012518x	*	*	20	392 x 306 x 375	4.00	0.045	100
4012718x	*	*	20	392 x 306 x 375	4.15	0.045	100
4018418x	*	*	20	392 x 306 x 375	4.02	0.045	20
4022718x	*	*	20	392 x 306 x 375	3.60	0.045	100
4022118x	*	*	20	392 x 306 x 375	3.50	0.045	100
4022218x	*	*	20	392 x 306 x 375	3.60	0.045	100
4022318x	*	*	20	392 x 306 x 375	4.00	0.045	100
4022418x	*	*	20	392 x 306 x 375	4.15	0.045	100
4022818x	*	*	20	392 x 306 x 375	4.02	0.045	100
4023218x	*	*	20	392 x 306 x 375	3.35	0.045	100
4023318x	*	*	20	392 x 306 x 375	3.40	0.045	100
4023418x	*	*	20	392 x 306 x 375	3.40	0.045	100
4023518x	*	*	20	392 x 306 x 375	4.00	0.045	100
4023618x	*	*	20	392 x 306 x 375	4.00	0.045	100
4023718x	*	*	20	392 x 306 x 375	4.15	0.045	100
40067069	100	580 x 380 x 176	200	591 x 391 x 385	8.10	0.089	200
40069069	100	580 x 380 x 176	200	591 x 391 x 385	14.00	0.089	200

* = Single packed - no shelf box

* Here noted as x, reference codes contain an additional digit (1-5) representing a language code.

Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Storage and handling:



With reservations for misprints and alterations in product specifications. Printed version not controlled.

Manufacturer: **ConvaTec Limited**
 First Avenue
 Address: Deeside Industrial Park
 Deeside, Flintshire, CH5 2NU, UK

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.
 The ConvaTec logo is a trademark of ConvaTec Inc.

