

UnoMeter™ Safeti™ Plus System

Hourly Diuresis Monitoring System (with Fixed Bag and KombiKon™ Sampling Port)



Intended use:

Urine meter with chamber for accurate measuring of hourly diuresis.
Fitted with bottom outlet. Designed to prevent infection.

List of product range:

REF	
158101310190	2.0L Fixed, straps and round hooks, 150 cm, free flow
158100510190	2.0L Fixed, straps and round hooks, 110 cm, free flow
158101210190	2.0L Fixed, straps, 150 cm, free flow
158100410190	2.0L Fixed, straps, 110 cm, free flow

Material:

Contains latex **No**
 Contains paper **Yes**
 Contains cotton **No**
 UnoMeter™ Safeti™ Plus . . . **Polyvinyl Chloride, ABS, Polypropylene, Polyethylene, Styrene Acrylonitrile**
 Airfilter **Teflon, non-woven Polyester backing**
 KombiKon™ **Polyvinyl Chloride, Polypropylene, Polyethylene, TPE**

Standardisation:

CE-marking **CE 0086**
 Quality system certified to . . **ISO 13485**
 GMDN **32072**

Sterilisation:

Sterile **Yes**
 Sterilised by **Ethylene Oxide**
 Shelf life **5 years from the date of sterilisation**

UnoMeter™ Safeti™ Plus System

Hourly Diuresis Monitoring System (with Fixed Bag and KombiKon™ Sampling Port)

UnoMeter™ Safeti™ Plus

Packaging:

REF	Shipping carton				MOQ pcs.
	Qty. Pcs.	Dimensions mm	Gross weight kg	Volume M ³	
158101310190	10	400x360x290	4.3	0.042	10
158100510190	10	400x360x290	4.3	0.042	10
158101210190	10	400x360x290	4.3	0.042	10
158100410190	10	400x360x290	4.3	0.042	10

Instruction for use:

Yes.

Storage and handling:

