





Aerosol Masks

Intended use:

Aerosol masks are intended for nebulization via a mask in order to administer medication via respiration. The product is intended for single use only.



Product range:

	REF*	Item	Size	Description
	93-100MM 93-100MMV (Vanilla scented mask)	Aerosol Mask	Adult	Adult Aerosol Mask (Elongated) with aerosol connector, Adjustable nose-clip and elastic strap
	93-220MM 93-220MMV (Vanilla scented mask)	Aerosol Mask	Pediatric	Pediatric Aerosol Mask (Elongated) with aerosol connector, Adjustable nose-clip and elastic strap
	93-260MM	Aerosol Mask	Infant	Infant Aerosol Mask with aerosol connector, Adjustable nose-clip and elastic strap
	93-110MM	Face Tent	Adult	Adult Face Tent with aerosol connector and elastic strap

* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Material specifications:

Contains latex? No
Contains phthalates? No
Contains bisphenols? No
Aerosol mask Polyvinyl chloride
Nose clip Aluminum
Elastic strap Polyester (Latex-free)
Aerosol connector High density polyethylene

Certifications:

CE-marking? Yes (CE 0120)
Product classification IIA
Quality system certified to EN ISO 13485:2012

Sterilization:

Sterile? No
Shelf life 5 years

Packaging specifications:

REF	Primary package	Secondary package/Polybag	Tertiary package/Shipper			Pallet
	Qty, pcs	Qty, pcs	Qty, pcs	Dimensions, mm (LxWxD)	Gross Wt., kg	Qty, pcs
93-100MM	1	10	50	400 x 300 x 240	1.91	1,600
93-100MMV	1	10	50	400 x 300 x 240	2.03	1,600
93-220MM	1	10	50	350 x 300 x 200	1.48	1,600
93-220MMV	1	10	50	350 x 300 x 200	1.28	1,600
93-260MM	1	10	50	350 x 300 x 200	1.40	1,600
93-110MM	1	10	50	350 x 300 x 200	1.48	900


EAN codes:

REF	Primary package	Secondary package	Tertiary package
93-100MM	(01)06945099302052	(01)26945099302056	(01)46945099302050
93-100MMV	(01)06945099302113	(01)26945099302117	(01)46945099302111
93-220MM	(01)06945099302069	(01)26945099302063	(01)46945099302067
93-220MMV	(01)06945099302120	(01)26945099302124	(01)46945099302128
93-260MM	(01)06945099302076	(01)26945099302070	(01)46945099302074
93-110MM	(01)06945099302106	(01)26945099302100	(01)46945099302104

Packaging material:

Primary package	Poly Bag - Low density polyethylene / Insert - 80g Paper
Secondary package	Poly Bag - Low density polyethylene / Sticker - Coated Paper
Tertiary package	Sticker - Coated Paper / Carton - Corrugated Board

Storage and handling:

	FRAGILE - Contents of the transport package are fragile and shall therefore be handled with care
	KEEP AWAY FROM RAIN - Transport package shall be kept away from rain
	THIS WAY UP - Indicates correct upright position of the transport package
	DO NOT USE HAND HOOKS - Hooks are prohibited for handling the transport package



Ningbo Shengyurui Medical Appliances Co., Ltd
 No. 138 Binhaisi Rd.
 Hangzhou Bay New Zone, 315336, Ningbo City
 Zhejiang Province, P.R. China
 Made in China

EC REP
Shanghai Intl. Holding Corp. GmbH
 Eiffestrasse 80, 20537 Hamburg
 Germany

Distributed by



ConvaTec Limited
 First Avenue, Deeside Industrial Park
 Deeside, Flintshire, CH5 2NU, UK
www.convatec.com

For further information on these products, please contact your local ConvaTec sales representative.

With reservations for misprints and alterations in product specifications. Printed version not controlled.
 All trade marks and logos are the property of a ConvaTec group company.

© 2017 ConvaTec Inc. AP-017302-MM, ver. FEB 2017