

# SimpaVac™

## Wound drainage system

**Intended use:**

SimpaVac™ is a low vacuum wound drainage system designed for post-operative wound drainage procedure. The product is intended for single use only.



**Product range:**

	REF*
SimpaVac™, 200 ml	4007118x

\* Here noted as x, reference codes contain an additional digit (1-5) representing a language code.  
 Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

**Material specifications:**

**Contains latex?** ..... No  
**Contains DEHP?** ..... No  
**O-ring** ..... Silicone  
**Container** ..... Polyethylene  
**Straps** ..... Polyethylene  
**Plug** ..... Polycarbonate

**Certifications:**

**CE-marking?** ..... Yes (CE 0086)  
**Product classification** ..... Class I (EU)  
**Quality system certified to** ..... ISO 13485  
**Global Medical Device Nomenclature (GMDN)** ..... 35824

**Sterilization:**

**Sterile?** ..... Yes  
**Sterilized by** ..... Ethylene oxide  
**Shelf life** ..... 5 years

**Packaging specifications:**

REF	Shipping carton				MOQ, pcs
	Qty, pcs	Dimensions, mm (LxWxD)	Gross Wt. kg	Volume m³	
4007118x	50	391 x 304 x 455	3.2	0.056	100

**Storage and handling:**




**ConvaTec Limited**  
 First Avenue  
 Deeside Industrial Park  
 Deeside, Flintshire, CH5 2NU, UK

Distributed by  
 **ConvaTec**  
 ConvaTec Limited  
 First Avenue, Deeside Industrial Park  
 Deeside, Flintshire, CH5 2NU, UK  
[www.convatec.com](http://www.convatec.com)

For further information on these products, please contact your local ConvaTec sales representative.

With reservations for misprints and alterations in product specifications. Printed version not controlled.  
 All trade marks and logos are the property of a ConvaTec group company.

© 2018 ConvaTec Inc. AP-018615-MM, ver. JAN 2018