




# Gas sampling lines

### Intended use:

Gas sampling lines are used to deliver gas samples to a monitoring device for the measurement of respiratory gases. The product is intended for single use only.



### List of product range:

Gas sampling line	REF*	Inner/outer diameter, mm	Length, m	Luer lock connectors
 895MM	895MM	1.68 / 3.33	3.0	Male/Male
	897MM	1.68 / 3.33	3.0	Male/Female
 897MM	892MM	1.27 / 3.33	3.0	Male/Male
	898MM	1.27 / 3.33	3.0	Male/Female

\* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

### Material specifications:

Contains latex? . . . . . No  
 Contains DEHP? . . . . . Yes  
 Gas sampling line . . . . . Polyvinyl chloride  
 Male/Female . . . . . Acrylic

### Certification:

CE-marking? . . . . . Yes (CE 0120)  
 Product classification . . . . . Class IIA (EU)  
 Quality system certified to . . . . . ISO 13485  
 Global Medical Device Nomenclature (GMDN) . . . . . 45566

### Sterilisation:

Sterile? . . . . . No  
 Shelf life . . . . . 5 years

### Packaging:

REF*	Shelf box		Shipping carton				MOQ, pcs.
	Qty., pcs.	Dimensions, mm	Qty., pcs.	Dimensions, mm	Gross Wt., kg	Volume, m <sup>3</sup>	
895MM	*	*	30	279 x 330 x 381	1.36	2.674	750
897MM	*	*	30	279 x 330 x 381	0.91	2.674	750
892MM	*	*	30	279 x 330 x 381	0.73	2.674	750
898MM	*	*	30	279 x 330 x 381	1.36	2.674	750

\* = Single packed - no shelf box

\* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

### Storage and handling:



*With reservations for misprints and alterations in product specifications. Printed version not controlled.*

Manufacturer: **Unomedical Inc.**  
 5701-1S, Ware Road  
 Address: McAllen, Texas 78503  
 USA

European Representative: **Unomedical a/s**  
 Address: Kongevejen 2  
 3460 Birkerød, Denmark

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.  
 The ConvaTec logo is a trademark of ConvaTec Inc.