


OP-Flex™ : Connecting Tubes


for OP-Flex™ suction handles

Intended use:

Connecting tubes for suction handles are intended to be used with a suction handle for the removal of blood and fluids during surgery. The product is intended for single use only.



List of product range:

Connecting tubes		REF**	Length, mm	CH/FG	Shore A
Two female connectors, one side cut to fit		1610018x	2,100	24	74
		1610218x	3,500	24	74

Material specifications:

Contains latex? No
Contains DEHP? Yes
Tube Polyvinyl chloride
Connectors Polyvinyl chloride

Certification:

CE-marking? Yes (CE 0086)
Product classification Class I (EU)
Quality system certified to ISO 13485
Global Medical Device Nomenclature (GMDN) 16779

Sterilisation:

Sterile? Yes
Sterilized by Ethylene oxide
Shelf life 5 years

Packaging:

REF**	Shipping carton*				MOQ, pcs.
	Qty, pcs.	Dimensions, mm	Gross Wt., kg	Volume, m ³	
1610018x	50	392 X 306 X 375	5.5	0.045	150
1610218x	50	392 X 306 X 375	6.8	0.045	300

* Single packed, no shelf box

** Here noted as x, reference codes contain an additional digit (1-5) representing a language code.
 Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Storage and handling:



With reservations for misprints and alterations in product specifications. Printed version not controlled.

Manufacturer: **ConvaTec Limited**
 First Avenue
 Address: Deeside Industrial Park
 Deeside, Flintshire, CH5 2NU, UK

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.
 The ConvaTec logo is a trademark of ConvaTec Inc.