


Siliconized Foley catheter

Tiemann/Coudé, 2-way, soft valve

Intended use:

The siliconized Foley catheter is a urinary catheter with a silicone fluid treated surface for easy passage through the urethral meatus.

The siliconized catheter is intended for short-term drainage of the bladder. The product is intended for single use only.



List of product range:

Siliconized Foley catheter, Tiemann/Coudé, 2-way, soft valve

CH size	5-10 ml/cc balloon	30 ml/cc balloon
12	MM51211210	MM51211230
14	MM51211410	MM51211430
16	MM51211610	MM51211630
18	MM51211810	MM51211830
20	MM51212010	MM51212030
22	MM51212210	MM51212230
24	MM51212410	MM51212430

Material specifications:

Contains latex? Yes
 Contains DEHP? No
 Catheter Latex
 Hard valve Rubber
 Sleeve Rubber
 Catheter surface Silicone emulsion

Certification:

CE-marking? Yes (CE 0123)
 Product classification Class IIB (EU)
 Quality system certified to ISO 13485
 Global Medical Device Nomenclature (GMDN). 34917

Sterilisation:

Sterile? Yes
 Sterilized by Gamma radiation
 Shelf life 5 years

Siliconized Foley catheter

Tiemann/Coudé, 2-way, soft valve

Packaging:

Siliconized Foley catheter, Tiemann/Coudé, 2-way, soft valve

REF*	Shelf box		Shipping carton				MOQ, pcs.
	Qty., pcs	Dimensions, mm	Qty., pcs	Dimensions, mm	Gross Wt., kg	Volume, m ³	
MM51211210	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM51211230	10	506 x 97 x 35	800	576 x 495 x 522	23	0.16	800
MM51211410	10	507 x 97 x 35	800	577 x 495 x 522	23	0.16	800
MM51211430	10	508 x 97 x 35	800	578 x 495 x 522	23	0.16	800
MM51211610	10	509 x 97 x 35	800	579 x 495 x 522	23	0.16	800
MM51211630	10	510 x 97 x 35	800	580 x 495 x 522	23	0.16	800
MM51211810	10	511 x 97 x 35	800	581 x 495 x 522	23	0.16	800
MM51211830	10	512 x 97 x 35	800	582 x 495 x 522	23	0.16	800
MM51212010	10	513 x 97 x 35	800	583 x 495 x 522	23	0.16	800
MM51212030	10	514 x 97 x 35	800	584 x 495 x 522	23	0.16	800
MM51212210	10	515 x 97 x 35	800	585 x 495 x 522	23	0.16	800
MM51212230	10	516 x 97 x 35	800	586 x 495 x 522	23	0.16	800
MM51212410	10	517 x 97 x 35	800	587 x 495 x 522	23	0.16	800
MM51212430	10	518 x 97 x 35	800	588 x 495 x 522	23	0.16	800

* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Storage and handling:



With reservations for misprints and alterations in product specifications. Printed version not controlled.

Manufacturer: **Unomedical Sdn Bhd**
(15290-P)
Address: Bakar Arang Industrial Estate
08000 Sungai Petani, Kedah, Malaysia

European Representative: **Unomedical a/s**
Address: Kongevejen 2
3460 Birkerød, Denmark

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.
The ConvaTec logo is a trademark of ConvaTec Inc.

