


Rectal Tubes

2 eyes

Intended use:

The rectal tube is inserted via the rectum to drain waste, perform lavage and/or administer medicine to lower intestine.
The product is intended for single use only.



List of product range:

	Effective length, mm	Tube Outer diameter of drain, mm									
		6.00	6.76	7.38	8.00	8.28	8.66	9.28	9.95	10.62	11.62
		CH18	CH20	CH22	CH24	CH25	CH26	CH28	CH30	CH32	CH35
		REF*	REF*	REF*	REF*	REF*	REF*	REF*	REF*	REF*	REF*
Rectal tubes 2 eyes	360	1400318x	1400618x	1400818x	1401018x	1401218x	1401418x	1401618x	1401818x	1402018x	1402218x

All REF: Shore A hardness: 80

* Here noted as x, reference codes contain an additional digit (1-5) representing a language code.
Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Material specifications:

Contains latex? No
Contains DEHP? No
Tubes Polyvinyl chloride
Connector Polyvinyl chloride

Certification:

CE-marking? Yes (CE 0086)
Product classification Class IIA (EU)
Quality system certified to ISO 13485
Global Medical Device Nomenclature (GMDN). 14227

Sterilisation:

Sterile? Yes
Sterilized by Ethylene oxide
Shelf life 5 years

Rectal Tubes

2 eyes

Packaging:

REF*	Shelf box		Shipping carton				MOQ, pcs.
	Qty., pcs	Dimensions, mm	Qty., pcs	Dimensions, mm	Gross Wt., kg	Volume m ³	
1400318x	100	145 x 93 x 446	800	391 x 304 x 455	9.3	0.054	1,600
1400618x	100	145 x 124 x 446	600	391 x 304 x 455	8.9	0.054	1,800
1400818x	100	145 x 124 x 446	600	391 x 304 x 455	9.7	0.054	1,800
1401018x	100	145 x 124 x 446	600	391 x 304 x 455	11.2	0.054	1,800
1401218x	100	145 x 124 x 446	600	391 x 304 x 455	11.2	0.054	1,800
1401418x	100	145 x 124 x 446	600	391 x 304 x 455	11.2	0.054	1,800
1401618x	50	145 x 93 x 446	400	391 x 304 x 455	8.5	0.054	1,200
1401818x	50	145 x 93 x 446	400	391 x 304 x 455	9.8	0.054	1,200
1402018x	50	145 x 93 x 446	400	391 x 304 x 455	10.6	0.054	4,000
1402218x	50	145 x 93 x 446	400	391 x 304 x 455	10.4	0.054	4,000

* Here noted as x, reference codes contain an additional digit (1-5) representing a language code.

Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Storage and handling:



With reservations for misprints and alterations in product specifications. Printed version not controlled.

Manufacturer: **ConvaTec Limited**
First Avenue
Address: Deeside Industrial Park
Deeside, Flintshire, CH5 2NU, UK

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.
The ConvaTec logo is a trademark of ConvaTec Inc.

