

A6 PreKon™

Closed urine drainage system with Silicone Foley catheter

Intended use:

Closed urine drainage system with preconnected Silicone Foley catheter for the collection of urine.

The product is intended for single use only.



List of product range:

A6 PreKon™ sterile closed urine drainage system with preconnected Silicone Foley catheter

REF*	Catheter size, CH	Item description
44 924 1S-xx	14	Sterile closed urine drainage system with KombiKon™ needle free sample port, 110 cm tubing, Drip chamber, Non-return valve (NRV), Urine drainage bag with 2000 ml bag capacity, Bottom Outlet/Cross valve Preconnected to an all Silicone Foley catheter, Standard, 2-way, hard valve
44 926 1S-xx	16	
44 928 1S-xx	18	

* Here noted as xx, reference codes contain two additional digits (10-50) representing a language code.

Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Material specifications:

Contains latex?	No
Contains DEHP?	Yes
Tube inlet/outlet, Drip chamber, Foil	Polyvinyl chloride
KombiKon™ needle-free sample port	Polyvinyl chloride, Polypropylene, High density polyethylene, Polyethylene, Poly(styrene-butadiene-styrene)
Cap	Low density polyethylene, High density polyethylene
Hanger, Tube clamp	Polypropylene, High density polyethylene
Front filter	Polyester
Bottom outlet	Polypropylene, High density polyethylene
Catheter tube and balloon	Silicone
Catheter hard valve	Polypropylene
Catheter sleeve	Polyvinyl chloride
Catheter tip	Silicone, Barium sulphate

Certification:

CE-marking?	Yes (0086)
Product classification	Class IIB (EU)
Quality system certified to	ISO 13485
Global Medical Device Nomenclature (GMDN)	58921

Sterilisation:

Sterile?	Yes
Sterilized by	Ethylene oxide
Shelf life	5 years

A6 PreKon™

Closed urine drainage system with Silicone Foley catheter

Packaging:

A6 PreKon™ sterile closed urine drainage system with preconnected Silicone Foley catheter

REF*	Shipping carton				MOQ, pcs.
	Qty., pcs.	Dimensions, mm	Gross Wt., kg	Volume, m ³	
44 924 1S-xx	20	582 x 250 x 234	4.1	0.034	840
44 926 1S-xx	20	582 x 250 x 234	4.1	0.034	840
44 928 1S-xx	20	582 x 250 x 234	4.1	0.034	840

* Here noted as xx, reference codes contain two additional digits (10-50) representing a language code.

Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Storage and handling:



With reservations for misprints and alterations in product specifications. Printed version not controlled.

Manufacturer: **ConvaTec Limited**
First Avenue
Address: Deeside Industrial Park
Deeside, Flintshire, CH5 2NU, UK

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.
The ConvaTec logo is a trademark of ConvaTec Inc.

