


# Drain-Fix™

Securement device

**Intended use:**

**Drain-Fix™** securement devices are intended to be used to secure drains, catheters and tubes ranging in sizes from CH5 to CH22. Drain-Fix™ securement devices are specially designed to stop migration, movement and accidental removal of the drain/catheter. The product is intended for single use only.



**List of product range:**

	REF*	Size, cm	Applicable drains, CH
Drain-Fix™, Small	680M	9.1 x 10.0	5-14
Drain-Fix™, Large	685M	9.1 x 10.0	12-22

### Material specifications:

<b>Contains latex?</b> .....	No
<b>Contains DEHP?</b> .....	No
<b>Carrier.</b> .....	Non-woven polypropylene
<b>Adhesive.</b> .....	Acrylic adhesive
<b>Adhesive on central part Drain-Fix™</b> ..	Hydrocolloid 0.6 mm
<b>Plastic component</b> .....	Thermoplastic elastomer
<b>Tapes</b> .....	Double sided adhesive tape
<b>Backing paper on lid</b> .....	Paper (siliconized on one side)
<b>Backing paper on skin contact area</b> .....	Paper SRP
<b>Pouch</b> .....	Gas-paper, Polyester/Polyethylene, Self-reinforcing polyethylene film

**Certification:**

CE-marking?	Yes (0086)
Product classification	Class I (EU)
Quality system certified to	ISO 13485
Global Medical Device Nomenclature (GMDN)	56630

### Sterilisation:

<b>Sterile?</b>	Yes
<b>Sterilized by</b>	Gamma radiation
<b>Shelf life</b>	3 years

## Packaging:

	REF*	Shelf box		Shipping carton			MOQ, pcs.
		Qty., pcs	Dimensions, mm	Qty., pcs	Dimensions, mm	Volume, m³	
Drain-Fix™, Small	680M	25	207 x 138 x 135	200	564 x 220 x 290	0.036	2,500
Drain-Fix™, Large	685M	20	207 x 138 x 135	160	564 x 220 x 290	0.036	2,000

\* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

### Storage and handling:



*With reservations for misprints and alterations in product specifications. Printed version not controlled.*

Manufacturer: **ConvaTec Limited**  
First Avenue  
Address: Deeside Industrial Park  
Deeside, Flintshire, CH5 2NU, UK

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.  
The ConvaTec logo is a trademark of ConvaTec Inc.