

Mini vacuum set

Wound drainage system

Intended use:

Low vacuum wound drainage systems for post operative removal of blood and body fluids. The product is intended for single use only.



List of product range:

	REF**	Ch size	Length of drain, mm	Inner/outer diameter of drain, mm	Length of perforation/ Stainless steel introducer, mm
Mini vacuum set with 30 ml bellows, Redon drain* with introducer and X-ray marking	1300118x	CH 06	500	1.1/2.0	70/55
	1300218x	CH 08		1.6/2.7	58/55

*Drains: Shore A 69

Material specifications:

Contains latex? No
Contains DEHP? Yes
Bellows Polyethylene
Drains Polyvinyl chloride

Certification:

CE-marking? Yes (CE 0086)
Product classification Class IIA (EU)
Quality system certified to ISO 13485
Global Medical Device Nomenclature (GMDN) 35824

Sterilisation:

Sterile? Yes
Sterilized by Ethylene oxide
Shelf life 5 years

Packaging:

REF**	Shelf box		Shipping carton				MOQ, pcs.
	Qty., pcs.	Dimensions, mm	Qty., pcs.	Dimensions, mm	Gross Wt., kg	Volume, m ³	
1300118x	50	297 x 190 x 366	200	391 x 304 x 743	4.90	0.092	400
1300218x	50	297 x 190 x 366	200	391 x 304 x 743	4.90	0.092	400

** Here noted as x, reference codes contain an additional digit (1-5) representing a language code.
 Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Storage and handling:



With reservations for misprints and alterations in product specifications. Printed version not controlled.

Manufacturer: **ConvaTec Limited**
 First Avenue
 Address: **Deeside Industrial Park**
 Deeside, Flintshire, CH5 2NU, UK

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.
 The ConvaTec logo is a trademark of ConvaTec Inc.