

Thyroidea set

Wound drainage system

Intended use:

Low-vacuum wound drainage system for post operative removal of blood and body fluids after thyroidectomy.
The product is intended for single use only.



List of product range:

Thyroidea set	REF*	CH
100 ml bellows providing initial vacuum of 140 mmHg /18 kPa 2 Ulmer drains: Shore A 69, Inner/outer diameter of drains 2.0/3.3 mm, Length of drains 500 mm Length of perforation 50 mm, 150 mm connecting tube with clamp Connector for optional use of 1 or 2 drains	1300418x	10

Material specifications:

Contains latex? No
Contains DEHP? Yes
Bellows Polyethylene
Drains and connecting tube Polyvinyl chloride
Trocar Stainless steel
Connector for 2 drains Silicone

Certification:

CE-marking? Yes (CE 0086)
Product classification Class IIA (EU)
Quality system certified to ISO 13485
Global Medical Device Nomenclature (GMDN) 35824

Sterilisation:

Sterile? Yes
Sterilized by Ethylene oxide
Shelf life 5 years

Packaging:

REF*	Shelf box		Shipping carton				MOQ. pcs.
	Qty., pcs	Dimensions. mm	Qty., pcs	Dimensions. mm	Gross Wt., kg	Volume m ³	
1300418x	10	294 x 125 x 360	60	391 x 304 x 743	5.50	0.092	180

* Here noted as x. reference codes contain an additional digit (1-5) representing a language code.
 Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Storage and handling:



Manufacturer: **ConvaTec Limited**
 First Avenue
 Address: Deeside Industrial Park
 Deeside, Flintshire, CH5 2NU, UK

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.
 The ConvaTec logo is a trademark of ConvaTec Inc.

With reservations for misprints and alterations in product specifications. Printed version not controlled.