

Product Data Sheet

# UnoMeter™ Abdo-Pressure™ IAP Monitoring System Kit/Stand-alone



**Intended use:**

Manometer for monitoring intra-abdominal pressure (IAP). Must be connected to a urethral Foley catheter (Kit/Stand-alone) and urine collection device (Stand-alone only).

**Product range:**

Description	REF
UnoMeter™ Abdo-Pressure™ IAP monitoring system kit (with UnoMeter™ Safeti™ Plus, 110cm tubing and round hook)	158100910190
UnoMeter™ Abdo-Pressure™ IAP monitoring system stand-alone	158100110190

**Material:**

Contains Latex? . . . . . No  
 Contains paper (packaging only)? . . . . . Yes  
 Contains cotton? . . . . . No  
 UnoMeter™ Abdo-Pressure™ urine meter. Polyethylene (Clamp), Polyvinyl Chloride (Tubing)  
 Airfilter . . . . . 1,2 micron Acrylic Copolymer hydrophobic membrane  
 UnoMeter™ Safeti™ Plus urine meter . . . . . Polyvinyl Chloride, ABS, Polypropylene,  
 Styrene Acrylonitrile Polytetrafluoroethylene (PTFE)  
 KombiKon™ needle-free sample port. Polyvinyl Chloride, Polypropylene, Polyethylene, TPE

**Standardisation:**

CE-marking . . . . . CE 0086  
 Quality system certified to . . . . . ISO 13485  
 GMDN . . . . . 47505  
 Medical Device Classification . . . . . Class IIA

**Sterilisation:**

Sterile . . . . . Yes  
 Sterilised by . . . . . Ethylene Oxide

**Packaging:**

REF	Shelf box		Shipping carton				MOQ pcs.
	Qty. pcs.	Dimensions mm	Qty. pcs.	Dimensions mm	Gross Wt. kg	Volume m <sup>3</sup>	
158100910190	na	na	12	460x380x290	5,6	0,051	120
158100110190	20	410x294x167	480	600x830x1020	70	0,53	20

**Instruction for use:**

No .....

**Storage and handling:**



TM indicates a trade mark of Unomedical a/s.  
The ConvaTec logo is a trademark of ConvaTec Inc.