

Siliconized Foley catheter

Standard, 3-way, soft valve/hard valve

Intended use:

The siliconized Foley catheter is a urinary catheter with a silicone fluid treated surface for easy passage through the urethral meatus. The siliconized Foley catheter is intended for short-, medium- and long-term drainage of the bladder. The product is intended for single use only.



List of product range:

Siliconized Foley catheter, Standard, 3-way, soft valve/hard valve

| CH size | 3-way SV 5-10 ml/cc balloon | 3-way SV 30 ml/cc balloon | 3-way HV 5-10 ml/cc balloon | 3-way HV 30 ml/cc balloon |
|---------|--------------------------------|------------------------------|--------------------------------|------------------------------|
| 16 | MM55111610 | MM55111630 | MM55121610 | MM55121630 |
| 18 | MM55111810 | MM55111830 | MM55121810 | MM55121830 |
| 20 | MM55112010 | MM55112030 | MM55122010 | MM55122030 |
| 22 | MM55112210 | MM55112230 | MM55122210 | MM55122230 |
| 24 | MM55112410 | MM55112430 | MM55122410 | MM55122430 |
| 26 | MM55112610 | MM55112630 | MM55122610 | MM55122630 |

Material specifications:

Contains latex? Yes
 Contains DEHP (Hard valve)? Yes
 Contains DEHP (Soft valve)? No
 Catheter Latex
 Hard valve Polypropylene
 Soft valve Rubber
 Sleeve (Hard valve) Polyvinyl chloride
 Sleeve (Soft valve) Rubber
 Catheter surface Silicone emulsion

Certification:

CE-marking? Yes (CE 0123)
 Product classification Class IIB (EU)
 Quality system certified to ISO 13485
 Global Medical Device Nomenclature (GMDN).34917

Sterilisation:

Sterile? Yes
 Sterilized by Gamma radiation
 Shelf life 5 years

Siliconized Foley catheter

Standard, 3-way, soft valve/hard valve

Packaging:

Siliconized Foley catheter, Standard, 3-way, soft valve/hard valve

| REF* | Shelf box | | Shipping carton | | | | MOQ, pcs. |
|------------|-----------|----------------|-----------------|-----------------|---------------|------------------------|-----------|
| | Qty., pcs | Dimensions, mm | Qty., pcs | Dimensions, mm | Gross Wt., kg | Volume, m ³ | |
| MM55111610 | 5 | 505 x 97 x 35 | 400 | 575 x 495 x 522 | 18.6 | 0.17 | 400 |
| MM55111630 | 5 | 505 x 97 x 35 | 400 | 575 x 495 x 522 | 18.6 | 0.17 | 400 |
| MM55111810 | 5 | 505 x 97 x 35 | 400 | 575 x 495 x 522 | 18.6 | 0.17 | 400 |
| MM55111830 | 5 | 505 x 97 x 35 | 400 | 575 x 495 x 522 | 18.6 | 0.17 | 400 |
| MM55112010 | 5 | 505 x 97 x 35 | 400 | 575 x 495 x 522 | 18.6 | 0.17 | 400 |
| MM55112030 | 5 | 505 x 97 x 35 | 400 | 575 x 495 x 522 | 18.6 | 0.17 | 400 |
| MM55112210 | 5 | 505 x 97 x 35 | 400 | 575 x 495 x 522 | 18.6 | 0.17 | 400 |
| MM55112230 | 5 | 505 x 97 x 35 | 400 | 575 x 495 x 522 | 18.6 | 0.17 | 400 |
| MM55112410 | 5 | 505 x 97 x 35 | 400 | 575 x 495 x 522 | 18.6 | 0.17 | 400 |
| MM55112430 | 5 | 505 x 97 x 35 | 400 | 575 x 495 x 522 | 18.6 | 0.17 | 400 |
| MM55112610 | 5 | 505 x 97 x 35 | 400 | 575 x 495 x 522 | 18.6 | 0.17 | 400 |
| MM55112630 | 5 | 505 x 97 x 35 | 400 | 575 x 495 x 522 | 18.6 | 0.17 | 400 |
| MM55121610 | 5 | 505 x 97 x 35 | 400 | 575 x 495 x 522 | 18.6 | 0.17 | 400 |
| MM55121630 | 5 | 505 x 97 x 35 | 400 | 575 x 495 x 522 | 18.6 | 0.17 | 400 |
| MM55121810 | 5 | 505 x 97 x 35 | 400 | 575 x 495 x 522 | 18.6 | 0.17 | 400 |
| MM55121830 | 5 | 505 x 97 x 35 | 400 | 575 x 495 x 522 | 18.6 | 0.17 | 400 |
| MM55122010 | 5 | 505 x 97 x 35 | 400 | 575 x 495 x 522 | 18.6 | 0.17 | 400 |
| MM55122030 | 5 | 505 x 97 x 35 | 400 | 575 x 495 x 522 | 18.6 | 0.17 | 400 |
| MM55122210 | 5 | 505 x 97 x 35 | 400 | 575 x 495 x 522 | 18.6 | 0.17 | 400 |
| MM55122230 | 5 | 505 x 97 x 35 | 400 | 575 x 495 x 522 | 18.6 | 0.17 | 400 |
| MM55122410 | 5 | 505 x 97 x 35 | 400 | 575 x 495 x 522 | 18.6 | 0.17 | 400 |
| MM55122430 | 5 | 505 x 97 x 35 | 400 | 575 x 495 x 522 | 18.6 | 0.17 | 400 |
| MM55122610 | 5 | 505 x 97 x 35 | 400 | 575 x 495 x 522 | 18.6 | 0.17 | 400 |
| MM55122630 | 5 | 505 x 97 x 35 | 400 | 575 x 495 x 522 | 18.6 | 0.17 | 400 |

* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Storage and handling:



With reservations for misprints and alterations in product specifications. Printed version not controlled.

Manufacturer: **Unomedical Sdn Bhd**
(15290-P)
Address: Bakar Arang Industrial Estate
08000 Sungai Petani, Kedah, Malaysia

European Representative: **Unomedical a/s**
Address: Kongevejen 2
3460 Birkerød, Denmark

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.
The ConvaTec logo is a trademark of ConvaTec Inc.

