

Careline™ leg bags

Sterile

Intended use:

Urine drainage bag with scale for collection of urine output. The device is not intended to have a measurement function. The scale is approximate and for indication purposes only. The product is intended for single use only.



List of product range:

Careline™ leg bags, Sterile

REF***	Volume ml	Tube length cm	KombiKon™ needle-free sample port/ Connector	Connector	Cut to fit	Cross valve	Lever tap	Textile film on backside of bag	Hospital Pack*	Community Pack**
12 618 1S-xx	500	8.5	Yes		No	Yes		No	Yes	
12 629 1S-xx	500	8.5	Yes		No		Yes	No	Yes	
12 623 1S-xx	500	29.0	Yes		No	Yes		No	Yes	
12 632 1S-xx	500	29.0	Yes		No		Yes	No	Yes	
08 613 1S-xx	500	8.5	Yes		No	Yes		Yes		Yes
08 620 1S-xx	500	29.0	Yes		No	Yes		Yes		Yes
12 613 1S-xx	500	8.5	Yes		No	Yes		Yes		Yes
12 612 1S-xx	500	8.5	Yes		No		Yes	Yes		Yes
12 620 1S-xx	500	29.0	Yes		No	Yes		Yes		Yes
12 621 1S-xx	500	29.0	Yes		No		Yes	Yes		Yes
24 111 1S-xx	750	10.0	No	Yes	No	Yes		Yes		Yes
24 151 1S-xx	750	50.0	No	Yes	Yes	Yes		Yes		Yes
24 610 1S-xx	750	8.5	Yes		No	Yes		Yes		Yes
24 612 1S-xx	750	9.0	Yes		No		Yes	Yes		Yes
24 620 1S-xx	750	29.0	Yes		No	Yes		Yes		Yes
24 621 1S-xx	750	29.0	Yes		No		Yes	Yes		Yes

*Includes 2 leg straps per leg bag, ** Includes 2 leg straps per shelf box (10 leg bags)

*** Here noted as xx, reference codes contain two additional digits (10-50) representing a language code. Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Material specifications:

- Contains latex?** No
- Contains DEHP?** Yes
- Contains cotton?** No
- Tube inlet/outlet, Foil.** Polyvinyl chloride
- Sample port/KombiKon™** . Polyvinyl chloride, Polypropylene, High density polyethylene, Polyethylene, Poly(styrene-butadiene-styrene)
- Regular connector.** Polyvinyl chloride
- Cross valve.** Polypropylene, High density polyethylene
- Cap** Low density polyethylene, High density polyethylene
- Tube on bottom outlet.** Silicone
- Lever Tap** Polyvinyl chloride, Acrylonitrile-butadiene-styrene, Silicone
- Films** Polyvinyl chloride
- Strap** Polyamide

Certification:

- CE-marking?** Yes (CE 0086)
- Product classification** Class I (EU)
- Quality system certified to** ISO 13485
- Global Medical Device Nomenclature (GMDN)** 58923

Sterilisation:

- Sterile?** Yes
- Sterilized by** Ethylene oxide
- Shelf life** 5 years

Careline™ leg bags

Sterile

Packaging:

REF*	Shelf box		Shipping carton				MOQ, pcs.
	Qty., pcs.	Dimensions, mm	Qty., pcs.	Dimensions, mm	Gross Wt., kg	Volume, m ³	
08 613 1S-xx	10	210 x 190 x 150	100	792 x 218 x 392	6.6	0.067	2,800
08 620 1S-xx	10	210 x 190 x 150	100	792 x 218 x 392	6.6	0.067	2,800
12 612 1S-xx	10	210 x 190 x 150	100	792 x 218 x 392	6.0	0.067	2,000
12 613 1S-xx	10	210 x 190 x 150	100	792 x 218 x 392	6.3	0.067	2,000
12 618 1S-xx	10	210 x 190 x 150	100	792 x 218 x 392	5.6	0.067	2,000
12 620 1S-xx	10	210 x 190 x 150	100	792 x 218 x 392	7.7	0.067	2,000
12 621 1S-xx	10	210 x 190 x 150	100	792 x 218 x 392	7.9	0.067	2,000
12 623 1S-xx	10	210 x 190 x 150	100	792 x 218 x 392	6.4	0.067	2,000
12 629 1S-xx	10	210 x 190 x 150	100	792 x 218 x 392	6.0	0.067	2,000
12 632 1S-xx	10	210 x 190 x 150	100	792 x 218 x 392	6.6	0.067	2,000
24 111 1S-xx	10	280 x 188 x 117	40	388 x 286 x 242	2.2	0.027	2,240
24 151 1S-xx	10	280 x 188 x 117	40	388 x 286 x 242	2.2	0.027	2,240
24 610 1S-xx	10	210 x 190 x 150	100	792 x 218 x 392	7.1	0.067	2,000
24 612 1S-xx	10	210 x 190 x 150	100	792 x 218 x 392	7.1	0.067	2,000
24 620 1S-xx	10	210 x 190 x 150	100	792 x 218 x 392	6.9	0.067	2,000
24 621 1S-xx	10	210 x 190 x 150	100	792 x 218 x 392	7.1	0.067	2,000

* Here noted as xx, reference codes contain two additional digits (10-50) representing a language code.

Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Storage and handling:



With reservations for misprints and alterations in product specifications. Printed version not controlled.

Manufacturer: **ConvaTec Limited**
 First Avenue
 Address: Deeside Industrial Park
 Deeside, Flintshire, CH5 2NU, UK

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.
 The ConvaTec logo is a trademark of ConvaTec Inc.

