


Tapered Yankauer

Suction handles and sets

Intended use:

The Tapered Yankauer is a suction handle intended for the removal of blood and fluids during surgery. The tapered Yankauer is designed with a narrower tip, used when space is limited and to increase the accuracy of suctioning. The product is intended for single use only.



List of product range:

	REF*	Connecting tube Shore A 74		Outer diameter of suction tip mm	Suction handle, length in mm	Vacuum control	Number of relieving eyes	CH/FG (handle)	
		Length in mm	CH/FG						
Tapered Yankauer, single bend	3402418x	2,100	24	3.8	210	No	0	12	Small
	3402518x	3,000	24	3.8	210	No	0	12	Small
	3302218x			3.8	210	No	0	12	Small

Material specifications:

Contains latex? No
Contains DEHP Yes
Handle Polyvinyl chloride
Tubing Polyvinyl chloride
Other components Polyvinyl chloride

Certification:

CE-marking? Yes (CE 0086)
Product classification Class IIA (EU)
Quality system certified to ISO 13485
Global Medical Device
Nomenclature (GMDN) 35917/36778

Sterilisation:

Sterile? Yes
Sterilized by Ethylene oxide
Shelf life 5 years

Packaging:

REF*	Shelf box		Shipping carton				MOQ, pcs.
	Qty., pcs	Dimensions, mm	Qty., pcs	Dimensions, mm	Gross Wt., kg	Volume, m ³	
3302218x	50	294 x 125 x 360	300	391 x 304 x 743	5.60	0.092	600
3402418x	*	*	50	392 x 306 x 340	6.00	0.043	150
3402518x	*	*	50	392 x 306 x 340	4.80	0.043	150

* = Single packed - no shelf box

* Here noted as x, reference codes contain an additional digit (1-5) representing a language code. Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Storage and handling:



With reservations for misprints and alterations in product specifications. Printed version not controlled.

Manufacturer: **ConvaTec Limited**
 First Avenue
 Address: Deeside Industrial Park
 Deeside, Flintshire, CH5 2NU, UK

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.
 The ConvaTec logo is a trademark of ConvaTec Inc.