

















Tracheal Tubes

Standard










Intended use:

Tracheal Tubes are inserted through the mouth or nose into the trachea and connected to ventilation equipment to secure an airway. The products are intended for single use only.


Product range:
Magill shaped Tracheal Tubes with HVLP* cuff

With Murphy eye REF	Without Murphy eye REF	With stylet and Murphy eye REF	Size		Length mm*	CRD mm**	Connector colour***
			I.D. mm	O.D. mm			
UM61110030	UM61120030		3.0	4.2	165	12	Red 
UM61110035	UM61120035		3.5	4.9	185	12	Red 
UM61110040	UM61120040		4.0	5.5	210	14	Green 
UM61110045	UM61120045		4.5	6.2	225	14	Green 
UM61110050	UM61120050	UM61110150	5.0	6.9	245	17	Blue 
UM61110055	UM61120055	UM61110155	5.5	7.5	275	17	Blue 
UM61110060	UM61120060	UM61110160	6.0	8.2	285	20	Blue 
UM61110065	UM61120065	UM61110165	6.5	8.8	295	20	Black 
UM61110070	UM61120070	UM61110170	7.0	9.6	310	25	Black 
UM61110075	UM61120075	UM61110175	7.5	10.2	315	25	White 
UM61110080	UM61120080	UM61110180	8.0	10.9	330	26	White 
UM61110085	UM61120085	UM61110185	8.5	11.5	330	26	Green 
UM61110090	UM61120090		9.0	12.1	330	28	Green 
UM61110095	UM61120095		9.5	12.7	330	28	Orange 
UM61110100	UM61120100		10.0	13.6	330	28	Orange 

Magill shaped Tracheal Tubes without HVLP* cuff

With Murphy eye REF	Without Murphy eye REF	Size		Length mm*	Connector colour***
		I.D. mm	O.D. mm		
UM61130020	UM61140020	2.0	2.9	140	Clear 
UM61130025	UM61140025	2.5	3.5	145	Clear 
UM61130030	UM61140030	3.0	4.2	165	Red 
UM61130035	UM61140035	3.5	4.9	185	Red 
UM61130040	UM61140040	4.0	5.5	210	Green 
UM61130045	UM61140045	4.5	6.2	225	Green 
UM61130050	UM61140050	5.0	6.9	245	Blue 
UM61130055	UM61140055	5.5	7.5	275	Blue 
UM61130060	UM61140060	6.0	8.2	285	Blue 

* HVLP = High Volume Low Pressure; Length is measured from tip to angle

** CRD = Cuff Resting Diameter

*** 15 mm ISO connector (ISO 8836)

Tracheal Tubes

Standard



Material specifications:

Contains latex? No
Contains phthalates? No
Tube Polyvinyl chloride
Cuff Polyvinyl chloride
X-ray line Polyvinyl chloride
Inflation line Polyvinyl chloride
Connector Polypropylene
Pilot balloon Polyvinyl chloride
Valve Polyvinyl chloride, Polypropylene,
 Acrylonitrile-butadiene-styrene, Silicone, Stainless steel
Stylet Aluminum
Sheath Polyvinyl chloride

Certifications:

CE-marking? Yes (0123)
Product classification Class IIA (EU)
Quality system certified to ISO 13485
Global Medical Device Nomenclature (GMDN) 46967

Sterilization:

Sterile? Yes
Sterilized by Ethylene oxide
Shelf life 5 years

Packaging specifications:

REF*	Sizes	Primary package/ Single peelpack		Secondary package/ Shelf box		Tertiary package/ Shipper				Pallet Qty pcs
		Qty pcs	Dimensions, mm (LxW)	Qty pcs	Dimensions, mm (LxWxD)	Qty pcs	Dimensions, mm (LxWxD)	Gross Wt. kg	Volume m ³	
Magill shaped Tracheal Tubes with HVLP cuff										
UM611100XX UM611200XX	3.0 to 4.5	1	271 x 90	10	280 x 125 x 90	100	470 x 295 x 270	5.0	0.037	1,800
UM611100XX UM611200XX	5.0 to 10.0	1	347 x 90	10	355 x 160 x 94	100	495 x 370 x 320	6.0	0.059	1,200
Magill shaped Tracheal Tubes without HVLP cuff										
UM611300XX UM611400XX	2.0 to 4.5	1	271 x 90	10	280 x 125 x 90	100	470 x 295 x 270	5.0	0.037	1,800
UM611300XX UM611400XX	5.0 to 6.0	1	347 x 90	10	355 x 160 x 94	100	495 x 370 x 320	6.0	0.059	1,600
Magill shaped Tracheal Tubes with preloaded stylet and HVLP cuff										
UM611101XX	5.0 to 8.5	1	347 x 90	10	355 x 160 x 94	100	495 x 370 x 320	7.0	0.059	1,200

* XX= size, e.g. 45 = 4.5 mm I.D.



Tracheal Tubes

Standard



EAN codes:

Magill shaped Tracheal Tubes with HVLP cuff and Murphy eye

REF	Primary package	Secondary package	Tertiary package
UM61110030	(01)30768455148173	(01)50768455148177	(01)70768455148171
UM61110035	(01)30768455148180	(01)50768455148184	(01)70768455148188
UM61110040	(01)30768455148197	(01)50768455148191	(01)70768455148195
UM61110045	(01)30768455148203	(01)50768455148207	(01)70768455148201
UM61110050	(01)30768455148210	(01)50768455148214	(01)70768455148218
UM61110055	(01)30768455148227	(01)50768455148221	(01)70768455148225
UM61110060	(01)30768455148234	(01)50768455148238	(01)70768455148232
UM61110065	(01)30768455148241	(01)50768455148245	(01)70768455148249
UM61110070	(01)30768455148258	(01)50768455148252	(01)70768455148256
UM61110075	(01)30768455148265	(01)50768455148269	(01)70768455148263
UM61110080	(01)30768455148272	(01)50768455148276	(01)70768455148270
UM61110085	(01)30768455148289	(01)50768455148283	(01)70768455148287
UM61110090	(01)30768455148296	(01)50768455148290	(01)70768455148294
UM61110095	(01)30768455148302	(01)50768455148306	(01)70768455148300
UM61110100	(01)30768455148319	(01)50768455148313	(01)70768455148317

Magill shaped Tracheal Tubes with HVLP cuff, without Murphy eye

REF	Primary package	Secondary package	Tertiary package
UM61120030	(01)30768455148326	(01)50768455148320	(01)70768455148324
UM61120035	(01)30768455148333	(01)50768455148337	(01)70768455148331
UM61120040	(01)30768455148340	(01)50768455148344	(01)70768455148348
UM61120045	(01)30768455148357	(01)50768455148351	(01)70768455148355
UM61120050	(01)30768455148364	(01)50768455148368	(01)70768455148362
UM61120055	(01)30768455148371	(01)50768455148375	(01)70768455148379
UM61120060	(01)30768455148388	(01)50768455148382	(01)70768455148386
UM61120065	(01)30768455148395	(01)50768455148399	(01)70768455148393
UM61120070	(01)30768455148401	(01)50768455148405	(01)70768455148409
UM61120075	(01)30768455148418	(01)50768455148412	(01)70768455148416
UM61120080	(01)30768455148425	(01)50768455148429	(01)70768455148423
UM61120085	(01)30768455148432	(01)50768455148436	(01)70768455148430
UM61120090	(01)30768455148449	(01)50768455148443	(01)70768455148447
UM61120095	(01)30768455148456	(01)50768455148450	(01)70768455148454
UM61120100	(01)30768455148463	(01)50768455148467	(01)70768455148461

Magill shaped Tracheal Tubes without HVLP cuff, with Murphy eye

REF	Primary package	Secondary package	Tertiary package
UM61130020	(01)30768455148500	(01)50768455148504	(01)70768455148508
UM61130025	(01)30768455148517	(01)50768455148511	(01)70768455148515
UM61130030	(01)30768455148524	(01)50768455148528	(01)70768455148522
UM61130035	(01)30768455148531	(01)50768455148535	(01)70768455148539
UM61130040	(01)30768455148548	(01)50768455148542	(01)70768455148546
UM61130045	(01)30768455148555	(01)50768455148559	(01)70768455148553
UM61130050	(01)30768455148562	(01)50768455148566	(01)70768455148560
UM61130055	(01)30768455148579	(01)50768455148573	(01)70768455148577
UM61130060	(01)30768455148586	(01)50768455148580	(01)70768455148584



Tracheal Tubes

Standard



EAN codes:

Magill shaped Tracheal Tubes without HVLP cuff and Murphy eye

REF	Primary package	Secondary package	Tertiary package
UM61140020	(01)30768455148593	(01)50768455148597	(01)70768455148591
UM61140025	(01)30768455148609	(01)50768455148603	(01)70768455148607
UM61140030	(01)30768455148616	(01)50768455148610	(01)70768455148614
UM61140035	(01)30768455148623	(01)50768455148627	(01)70768455148621
UM61140040	(01)30768455148630	(01)50768455148634	(01)70768455148638
UM61140045	(01)30768455148647	(01)50768455148641	(01)70768455148645
UM61140050	(01)30768455148654	(01)50768455148658	(01)70768455148652
UM61140055	(01)30768455148661	(01)50768455148665	(01)70768455148669
UM61140060	(01)30768455148678	(01)50768455148672	(01)70768455148676





Magill shaped Tracheal Tubes with preloaded stylet and HVLP cuff

REF	Primary package	Secondary package	Tertiary package
UM61110150	(01)30768455154853	(01)50768455154857	(01)70768455154851
UM61110155	(01)30768455154860	(01)50768455154864	(01)70768455154868
UM61110160	(01)30768455154877	(01)50768455154871	(01)70768455154875
UM61110165	(01)30768455154884	(01)50768455154888	(01)70768455154882
UM61110170	(01)30768455154891	(01)50768455154895	(01)70768455154899
UM61110175	(01)30768455154907	(01)50768455154901	(01)70768455154905
UM61110180	(01)30768455154914	(01)50768455154918	(01)70768455154912
UM61110185	(01)30768455154921	(01)50768455154925	(01)70768455154929

Packaging material:

Primary package	Peelback pouch - Uncoated paper, XPP film
Secondary package	Ivory cardboard
Tertiary package	Corrugated cardboard

Storage and handling:

	FRAGILE - Contents of the transport package are fragile and shall therefore be handled with care
	KEEP DRY
	THIS WAY UP - Indicates correct upright position of the transport package
	DO NOT USE HAND HOOKS - Hooks are prohibited for handling the transport package


Well Lead Medical Co.,Ltd.
 C-4 Jinhu Industrial Estate
 Hualong, Panyu, Guangzhou
 511434 P.R.China

Made in China



Shanghai Intl. Holding Corp. GmbH
 Eiffestrasse 80, 20537 Hamburg
 Germany

Distributed by



ConvaTec Limited
 First Avenue, Deeside Industrial Park
 Deeside, Flintshire, CH5 2NU, UK
www.convatec.com

For further information on these products, please contact your local ConvaTec sales representative.

With reservations for misprints and alterations in product specifications. Printed version not controlled.

All trade marks and logos are the property of a ConvaTec group company.

© 2018 ConvaTec Inc. AP-019246-MM, ver. JUL 2018