








Getting the benefits of me+™ is easy.

Simply call 0800 467 866 or visit www.convatec.co.uk/meplus

More than just great products — me+ brings you the tools and advice to help you make life with a stoma completely your own.

	Code	Description / Dimensions	Qty
	422160	ease™ Thin strips 3cm x 11 cm	60 pcs
	422161	ease™ Thin strips 5cm x 11 cm	30 pcs
	422162	ease™ Thin strips 5cm x 22 cm	20 pcs
	422163 (US and Brazil)	ease™ Thin curve 3cm x 9 cm	40 pcs
	422164	ease™ Thin curve 3cm x 9 cm	60 pcs

Support and advice every step of the way

To find out more about our ostomy solutions or if you'd like to talk to us, please get in touch on

XXXX XXX XXX

or visit www.convatec.com/ostomy

© 2017 ConvaTec Inc.

™ indicates a trademark of ConvaTec group company. AP-018293-MM



Distributed by



ease™

Ostomy care,
as individual
as you are

Lucy

Enjoying her weekly
coffee catch up

*Model portrayal

Strips

Thin flexible strips designed to adapt
to your body and move with you.

ease™ strips overlap the edge of the skin
barrier for extra security.

 ConvaTec

Flexible

Skin
friendly

Easy to
apply

ease™

Strips

ConvaTec

Strips

Helps give peace of mind to people living with an ostomy.

Thin flexible strips designed to adapt to your body and move with you. ease™ strips overlap the edge of the skin barrier for extra security.



Security


Adhesive backing keeps the strip in place and provides security.

If the edge of the barrier expands and becomes sticky, the strips cover this to prevent it from sticking to clothing.



Flexible and thin


Thin (0.44 mm) so it moves with the body.



Easy to apply and remove

Split paper backing allows half of the backing to be removed and applied to the edge of the barrier.

The remaining backing can then be removed once the strip is in situ.





Skin friendly

Made with skin friendly hydrocolloid.

Latex free.

Covers the edge of the barrier, so if this dries out and becomes hard it can help to reduce skin irritation.




Water resistant, yet breathable

Water resistant top layer so the strips can be used when bathing and swimming.

Latex Free

**Breathable
polyurethane
outer**

**Absorbs
Moisture**



Range of sizes

A range of sizes available to meet different individual's needs.



3 x 9 cm



3 x 11 cm



5 x 11 cm



5 x 22 cm

Please see back page for order codes