

# Guedel Airways

**Intended use:**

Guedel airways are intended to maintain the patency of the respiratory passages through the oral cavity and pharynx after general anaesthesia and in patients unconscious for other reasons.  
The product is intended for single use only.



**List of product range:**

REF*	Description	Size	Length, mm	Color coding
2203718x	Guedel Airways size 000	000	40	Pink
2202418x	Guedel Airways size 00	00	50	Blue
2202318x	Guedel Airways size 0	0	60	Black
2201918x	Guedel Airways size 1	1	70	White
2202018x	Guedel Airways size 2	2	80	Green
2202118x	Guedel Airways size 3	3	90	Yellow
2202218x	Guedel Airways size 4	4	100	Red
2205418x	Guedel Airways size 5	5	110	Light blue
2205518x	Guedel Airways size 6	6	120	Orange

\* Here noted as x, reference codes contain an additional digit (1-5) representing a language code.  
Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

**Material specifications:**

**Contains latex?** . . . . . No  
**Contains DEHP** . . . . . No  
**Main tube** . . . . . Low density polyethylene  
**Bite Reinforcement** . . . . . Polyoxymethylene

**Certification:**

**CE-marking?** . . . . . Yes (CE 0086)  
**Product classification** . . . . . Class I (EU)  
**Quality system certified to** . . . . . ISO 13485  
**Global Medical Device Nomenclature (GMDN)**. . . . .42424

**Sterilisation:**

**Sterile?** . . . . . Yes  
**Sterilized by** . . . . . Ethylene oxide  
**Shelf life** . . . . . 5 years

# Guedel Airways

## Packaging:

REF*	Shelf box		Shipping carton				MOQ, pcs.
	Qty., pcs	Dimensions, mm	Qty., pcs	Dimensions, mm	Gross Wt., kg	Volume, m <sup>3</sup>	
2201918x	100	294 x 125 x 360	300	392 x 306 x 375	7.9	0.045	600
2202018x	50	294 x 125 x 360	300	391 x 304 x 743	7.9	0.092	300
2202118x	50	294 x 125 x 360	300	391 x 304 x 743	8.5	0.092	600
2202218x	50	297 x 190 x 366	200	391 x 304 x 743	6.2	0.092	600
2202318x	100	294 x 125 x 360	300	392 x 306 x 375	6.9	0.045	900
2202418x	100	294 x 125 x 360	300	392 x 306 x 375	6.9	0.045	900
2203718x	100	294 x 125 x 360	300	392 x 306 x 375	2.0	0.045	900
2205418x	50	297 x 190 x 366	200	391 x 304 x 743	6.2	0.092	200
2205518x	50	297 x 190 x 366	200	391 x 304 x 743	6.2	0.092	200

\* Here noted as x, reference codes contain an additional digit (1-5) representing a language code.

Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

## Storage and handling:



*With reservations for misprints and alterations in product specifications. Printed version not controlled.*

Manufacturer: **ConvaTec Limited**  
First Avenue  
Address: Deeside Industrial Park  
Deeside, Flintshire, CH5 2NU, UK

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.  
The ConvaTec logo is a trademark of ConvaTec Inc.

