


Lubricious coated Foley catheter

Standard, 2-way and 3-way, hard valve

Intended use:

The lubricious coated Foley catheter is a urinary catheter intended for medium-term drainage of the bladder. The product is intended for single use only



List of product range:

Lubricious coated Foley catheter, Standard, 2-way and 3-way, hard valve

	Pediatric	Standard, 2-way	Standard, 2-way	Standard, 3-way	Standard, 3-way
CH size	3-5 ml/cc	5-10 ml/cc	30 ml/cc	5-10 ml/cc	30 ml/cc
6	MM53180605				
8	MM53180805				
10	MM53181005				
12		MM51181210	MM51181230		
14		MM51181410	MM51181430		
16		MM51181610	MM51181630	MM55181610	MM55181630
18		MM51181810	MM51181830	MM55181810	MM55181830
20		MM51182010	MM51182030	MM55182010	MM55182030
22		MM51182210	MM51182230	MM55182210	MM55182230
24		MM51182410	MM51182430	MM55182410	MM55182430
26		MM51182610	MM51182630	MM55182610	MM55182630
28		MM51182810	MM51182830		
30		MM51183010	MM51183030		

* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Material specifications:

Contains latex? Yes
 Contains DEHP? Yes
 Catheter tube and balloon Latex
 Hard valve Polypropylene
 Sleeve Polyvinyl chloride
 Catheter surface Hydrophilic coating

Certification:

CE-marking? Yes (CE 0123)
 Product classification Class IIB (EU)
 Quality system certified to ISO 13485
 Global Medical Device
 Nomenclature (GMDN) 34917 (2-way), 32331 (3-way)

Sterilisation:

Sterile? Yes
 Sterilized by Gamma radiation
 Shelf life 5 years

Lubricious coated Foley catheter

Standard, 2-way and 3-way, hard valve

Packaging:

Lubricious coated Foley catheter, Standard, 2-way and 3-way, hard valve

REF*	Shelf box		Shipping carton				MOQ, pcs.
	Qty., pcs	Dimensions, mm	Qty., pcs	Dimensions, mm	Gross Wt., kg	Volume, m ³	
MM53180605	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM53180805	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM53181005	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM51181210	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM51181410	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM51181610	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM51181810	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM51182010	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM51182210	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM51182410	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM51182610	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM51182810	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM51183010	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM51181230	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM51181430	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM51181630	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM51181830	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM51182030	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM51182230	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM51182430	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM51182630	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM51182830	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM51183030	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM55181610	5	505 x 97 x 35	400	575 x 495 x 522	23	0.16	400
MM55181810	5	505 x 97 x 35	400	575 x 495 x 522	23	0.16	800
MM55182010	5	505 x 97 x 35	400	575 x 495 x 522	23	0.16	800
MM55182210	5	505 x 97 x 35	400	575 x 495 x 522	23	0.16	800
MM55182410	5	505 x 97 x 35	400	575 x 495 x 522	23	0.16	800
MM55182610	5	505 x 97 x 35	400	575 x 495 x 522	23	0.16	800
MM55181630	5	505 x 97 x 35	400	575 x 495 x 522	23	0.16	800
MM55181830	5	505 x 97 x 35	400	575 x 495 x 522	23	0.16	800
MM55182030	5	505 x 97 x 35	400	575 x 495 x 522	23	0.16	800
MM55182230	5	505 x 97 x 35	400	575 x 495 x 522	23	0.16	800
MM55182430	5	505 x 97 x 35	400	575 x 495 x 522	23	0.16	800
MM55182630	5	505 x 97 x 35	400	575 x 495 x 522	23	0.16	800

* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Storage and handling:



With reservations for misprints and alterations in product specifications. Printed version not controlled.

Manufacturer: **Unomedical Sdn Bhd**
(15290-P)
Address: Bakar Arang Industrial Estate
08000 Sungai Petani, Kedah, Malaysia

European Representative: **Unomedical a/s**
Address: Kongevejen 2
3460 Birkerød, Denmark

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.
The ConvaTec logo is a trademark of ConvaTec Inc.

