

## 3-in-1 Emergency Masks

**Intended use:**


**Emergency masks** can be used for delivery of oxygen in emergency medicine. 3-in-1 emergency masks contain components to be used as high concentration non-rebreather, high concentration rebreather or medium concentration mask.

The product is intended for single use only.



**List of product range:**

### 3-in-1 Emergency masks, Elongated (under the chin)

	REF*	Size	Description	m	Type	Connector
	182MM	Adult	3-in-1 Emergency mask with adjustable side holes, Reservoir bag, with tubing	2.1	Standard lumen	Standard
	282MM	Pediatric	3-in-1 Emergency mask with adjustable side holes, Reservoir bag, with tubing	2.1	Standard lumen	Standard

### Material specifications:

<b>Contains latex?</b>	No
<b>Contains DEHP?</b>	Yes
<b>Mask</b>	Polyvinyl chloride
<b>Nose clip</b>	Aluminum
<b>Elastic strap</b>	(latex-free) Polyester
<b>Tubing</b>	Polyvinyl chloride
<b>Connectors</b>	(Connector) High density polyethylene (Bushing) Polyvinyl chloride
<b>Reservoir bag</b>	Polyvinyl chloride

## Certification:

<b>CE-marking?</b>	Yes (CE 0120)
<b>Product classification</b>	Class IIA (EU)
<b>Quality system certified to</b>	ISO 13485
<b>Global Medical Device Nomenclature (GMDN)</b>	12855

## Sterilisation:

**Sterile?** ..... No  
**Shelf life** ..... 5 years

## Packaging:

REF*	Shelf box		Shipping carton				MOQ, pcs.
	Qty., pcs	Dimensions, mm	Qty., pcs	Dimensions, mm	Gross Wt., kg	Volume, m³	
182MM	*	*	50	356 x 406 x 305	4.94	0.044	1,250
282MM	*	*	50	356 x 406 x 305	4.60	0.044	1,250

\* = Single packed - no shelf box

\* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

### Storage and handling:



*With reservations for misprints and alterations in product specifications. Printed version not controlled.*

Manufacturer: **Unomedical Inc.**  
5701-1S, Ware Road  
Address: McAllen, Texas 78503  
USA

European Representative: **Unomedical a/s**

Address: Kongevejen 2  
3460 Birkerød, Denmark

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.  
The ConvaTec logo is a trademark of ConvaTec Inc.