

A1

Basic urine drainage bags

Intended use:

Urine drainage bag for the collection and monitoring of urine output.
The product is intended for single use only.



List of product range:

REF*	Capacity, ml	Tube length, cm
33 001 0-xx	1,500	10
33 003 0-xx	1,500	35
33 005 0-xx	1,500	90
35 105 0-xx	2,000	90
35 030 3-xx	2,000	90

Material specifications:

Contains latex? No
 Contains DEHP? Yes
 Film Polyvinyl chloride
 Tube Polyvinyl chloride
 Connector Polyvinyl chloride
 Cap Low density polyethylene, High density polyethylene

Certification:

CE-marking? Yes
 Product classification Class I (EU)
 Quality system certified to ISO 13485
 Global Medical Device Nomenclature (GMDN) 58918

Sterilisation:

Sterile? No
 Shelf life 5 years

Packaging:

REF*	Shipping carton				MOQ, pcs.
	Qty., pcs.	Dimensions, mm	Gross Wt., kg	Volume, m ³	
33 001 0-xx	500	375 x 280 x 278	9.00	0.029	28,000
33 003 0-xx	250	375 x 280 x 180	5.80	0.019	6,000
33 005 0-xx	250	388 x 286 x 242	7.95	0.027	12,000
35 105 0-xx	250	375 x 280 x 278	9.80	0.029	14,000
35 030 3-xx	250	375 x 280 x 242	8.70	0.027	14,000

* Here noted as xx, reference codes contain two additional digits (10-50) representing a language code.

Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Storage and handling:



With reservations for misprints and alterations in product specifications. Printed version not controlled.

Manufacturer: **ConvaTec Limited**
 First Avenue
 Address: Deeside Industrial Park
 Deeside, Flintshire, CH5 2NU, UK

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.
 The ConvaTec logo is a trademark of ConvaTec Inc.