


# Flattened Yankauer

## Suction handles and sets

### Intended use:

The flattened Yankauer is a suction handle intended for the removal of blood and fluids during surgery. The flattened Yankauer is designed with a narrower tip, used when space is limited and to increase the accuracy of suctioning. The product is intended for single use only.



### List of product range:

	REF*	Connecting tube Shore A 74		Inner/outer diameter of suction, mm	Suction handle, length in mm	Vacuum control	Number of relieving eyes	CH/FG (handle)	
		Length in mm	CH/FG						
Flattened Yankauer, open tip, single bend	3301418x			4.6/6.5	210	No	0	20	Medium
	3400718x	2,100	24	4.6/6.5	210	No	0	20	Medium
	3400818x	3,000	24	4.6/6.5	210	No	0	20	Medium
	3402918x	2,100	30	8.0/10.0	210	No	0	30	X-Large

### Material specifications:

Contains latex? ..... No  
 Contains DEHP ..... Yes  
 Handle ..... Polyvinyl chloride  
 Tubing ..... Polyvinyl chloride  
 Other components ..... Polyvinyl chloride

### Certification:

CE-marking? ..... Yes (CE 0086)  
 Product classification ..... Class IIA (EU)  
 Quality system certified to ..... ISO 13485  
 Global Medical Device  
 Nomenclature (GMDN) ..... 35917/36778

### Sterilisation:

Sterile? ..... Yes  
 Sterilized by ..... Ethylene oxide  
 Shelf life ..... 5 years

### Packaging:

REF*	Shelf box		Shipping carton				MOQ, pcs.
	Qty., pcs	Dimensions, mm	Qty., pcs	Dimensions, mm	Gross Wt., kg	Volume, m <sup>3</sup>	
3301418x	50	294 x 125 x 360	300	391 x 304 x 743	3.90	0.049	600
3400718x	*	*	50	392 x 306 x 340	6.00	0.043	600
3400818x	*	*	50	392 x 306 x 340	6.40	0.049	150
3402918x	*	*	50	392 x 306 x 375	6.70	0.045	150

\* = Single packed - no shelf box

\* Here noted as x, reference codes contain an additional digit (1-5) representing a language code.

Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

### Storage and handling:



*With reservations for misprints and alterations in product specifications. Printed version not controlled.*

Manufacturer: **ConvaTec Limited**  
 First Avenue  
 Address: Deeside Industrial Park  
 Deeside, Flintshire, CH5 2NU, UK

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.  
 The ConvaTec logo is a trademark of ConvaTec Inc.