


# insuflon™

### Intended use:

The product is indicated for subcutaneous infusion of medication administered by injection.



### Product range:

Product code	Description	Cannula length	EAN	Pcs per inner box	Pcs per outer box
insuflon™					
98-000-2612	insuflon™ port	18 mm	05705244001756	10	200

### Product specifications:

**Insertion angle** . . . . . 20-45 degrees  
**Soft cannula length** . . . . . 18 mm  
**Soft cannula outer diameter** . . . . . ~ 0,68 mm  
**Introducer needle** . . . . . 27 gauge

### Material specifications:

**PVC free** . . . . . Yes  
**Made with natural rubber latex** . . . . . No  
**Raw materials or derivation of animal origin** . . . . . No  
**The product contains phthalates** . . . . . No  
**Introducer & connector needles may contain** . . ≤ 10% nickel

### Certifications:

**CE-marking** . . . . . Yes  
**Product classification** . . . . . Class IIB (EU)  
**Quality system certified to** . . . . . ISO 13485  
**Tariff code** . . . . . HS9018390000  
**GMDN** . . . . . 17180






### Sterilization:

**Sterile** . . . . . Yes  
**Sterilized by** . . . . . Ethylene oxide  
**Shelf life** . . . . . 5 years  
**Usage time** . . . . . Maximum 72 hours

### Drug compatibility/stability tested with:

Insulin

### Storage and handling:

	<b>AVOID DIRECT SUNLIGHT</b> - Keep the product away from direct sunlight and any humid environment
	<b>KEEP DRY</b> - Store in a dry place at room temperature
	<b>DO NOT RE-USE</b>
	<b>DO NOT USE</b> - If package is damaged
	<b>CAUTION</b> - See Instruction for Use