


Flexi-Seal[®] Faecal Collector

Collection Means Protection

For non-ambulatory patients who are incontinent of loose or semi-formed stool.



Flexi-Seal[®] Faecal Collector features a skin-friendly hydrocolloid adhesive with a tape border, plus a tapered pouch with a drainage port and deodorising filter.

Problems With Pads

- No odour containment
- Frequent changing, moving, and disturbing patients
- Need skin care products to protect skin
- Risk of skin erythema, maceration, and breakdown
- Extensive cleaning of patients

Advantages of Faecal Collection


- Helps protect skin integrity—Averts breakdown/maceration of perianal skin
- Helps minimise exposure to faeces
 - Facilitates measurement of output
 - Maintains patient comfort and care
 - May remain in place up to 48 hours
 - Can hold up to 1 Litre of liquid/semi-formed stool

Flexi-Seal[®] Faecal Collector

Change the Way You Change

To benefit the patient, the staff, and the hospital environment

- Easy to prepare — apply with minimal equipment: soft cloth for cleansing, skin cleanser, gloves, and scissors
- Designed to offer security and skin protection
- Equipped with a tap for use with a bedside drain unit for liquid stool, or a tail closure for semi-formed stool
- Does not contain latex



Product Ordering

		Code	NHS Code	
When Flexi-Seal [®] Faecal Management Systems [™] are not an option	Flexi-Seal[®] Faecal Collector With Hydrocolloid Adhesive	650078	GCT077	10 Collectors per box
For bedridden patients with liquid or semi-liquid stool	Flexi-Seal[®] CONTROL[™] Faecal Management System	411107	FWN101	

To learn more call

0800 289 738 (UK) or 1800 946 938 (ROI)

Mon - Fri, 9:00 AM - 5:00 PM

www.convatec.com



Please see Flexi-Seal[®] Faecal Collector With Hydrocolloid Adhesive package insert for complete instructions for use.

Please see Flexi-Seal[®] CONTROL[™] package insert for complete instructions for use.