


All Silicone Foley catheter

Standard, 3-way, hard valve

Intended use:

The all Silicone Foley catheter is a urinary catheter intended for short as well as long term drainage of the bladder. For routine and post-operative drainage. The product is intended for single use only.



List of product range:

All Silicone Foley catheter, Standard, 3-way, hard valve

CH size	5-10 ml/cc balloon	30 ml/cc balloon
18	MM45151805	MM45151830
20	MM45152005	MM45152030
22	MM45152205	MM45152230
24	MM45152405	MM45152430

Material specifications:

Contains latex? No
 Contains DEHP? Yes
 Catheter tube and balloon Silicone
 Valve Polypropylene
 Tip Silicone, Barium sulphate
 Sleeve Polyvinyl chloride

Certification:

CE-marking? Yes (CE 0123)
 Product classification Class IIB (EU)
 Quality system certified to ISO 13485
 Global Medical Device Nomenclature (GMDN). 32331

Sterilisation:

Sterile? Yes
 Sterilized by Ethylene oxide
 Shelf life 5 years

All Silicone Foley catheter

Standard, 3-way, hard valve

Packaging:

All Silicone Foley catheter, Standard, 3-way, hard valve

REF*	Shelf box		Shipping carton				MOQ, pcs.
	Qty., pcs	Dimensions, mm	Qty., pcs	Dimensions, mm	Gross Wt., kg	Volume, m ³	
MM45151805	10	500 x 92 x 52	500	580 x 500 x 527	19.00	0.16	500
MM45151830	10	500 x 92 x 52	500	580 x 500 x 527	19.00	0.16	500
MM45152005	10	500 x 92 x 52	500	580 x 500 x 527	23.00	0.16	500
MM45152030	10	500 x 92 x 52	500	580 x 500 x 527	23.05	0.16	500
MM45152205	10	500 x 92 x 52	500	580 x 500 x 527	24.00	0.16	500
MM45152230	10	500 x 92 x 52	500	580 x 500 x 527	24.00	0.16	500
MM45152405	10	500 x 92 x 52	500	580 x 500 x 527	25.09	0.16	500
MM45152430	10	500 x 92 x 52	500	580 x 500 x 527	25.45	0.16	500

* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Storage and handling:



With reservations for misprints and alterations in product specifications. Printed version not controlled.

Manufacturer: **Unomedical Sdn Bhd**
(15290-P)
Address: Bakar Arang Industrial Estate
08000 Sungai Petani, Kedah, Malaysia

European Representative: **Unomedical a/s**
Address: Kongevejen 2
3460 Birkerød, Denmark

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.
The ConvaTec logo is a trademark of ConvaTec Inc.

