

A6

Sterile closed system bags

Intended use:

The A6 sterile closed system bag is intended for collection of urine. The product is intended for single use only.



List of product range:

Sterile closed system bag with drip chamber and non-return valve (NRV), 2000 ml

REF*	Tube length cm	Standard sample port (for use with needle)	KombiKon™ needle-free sample port	Tube clamp	Drip chamber	Capacity ml	Integrated hanger	Bottom outlet/ Cross valve
41 049 1S-xx	110	Yes	-	Yes	Yes	2,000	Yes	Yes
41 149 1S-xx	110	-	Yes	Yes	Yes	2,000	Yes	Yes

Material specifications:

Contains latex? No
Contains DEHP? Yes
Tube inlet/outlet, Drip chamber, Foil Polyvinyl chloride
Sample port/KombiKon™ . Polyvinyl chloride, Polypropylene, High density polyethylene, Polyethylene, Poly(styrene-butadiene-styrene)
Cap Low density polyethylene, High density polyethylene
Hanger, Tube clamp Polypropylene, High density polyethylene
Front filter. Polyester
Bottom outlet Polypropylene, High density polyethylene

Certification:

CE-marking? Yes (0086)
Product classification Class I (EU)
Quality system certified to ISO 13485
Global Medical Device Nomenclature (GMDN) 58921

Sterilisation:

Sterile? Yes
Sterilized by Ethylene oxide
Shelf life 5 years

Packaging:

REF*	Shipping carton				MOQ, pcs.
	Qty., pcs.	Dimensions, mm	Gross Wt., kg	Volume, m ³	
41 049 1S-xx	20	375 x 280 x 278	3.65	0.029	960
41 149 1S-xx	20	375 x 280 x 278	4.02	0.029	960

* Here noted as xx, reference codes contain two additional digits (10-50) representing a language code. Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Storage and handling:



With reservations for misprints and alterations in product specifications. Printed version not controlled.

Manufacturer: **ConvaTec Limited**
 First Avenue
 Address: Deeside Industrial Park
 Deeside, Flintshire, CH5 2NU, UK

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.
 The ConvaTec logo is a trademark of ConvaTec Inc.