


# Siliconized Foley catheter

## Tiemann/Coudé, 2-way, soft valve

**Intended use:**

The siliconized Foley catheter is a urinary catheter with a silicone fluid treated surface for easy passage through the urethral meatus. The siliconized catheter is intended for short-term drainage of the bladder. The product is intended for single use only.



**List of product range:**

Siliconized Foley catheter, Tiemann/Coudé, 2-way, soft valve

CH size	5-10 ml/cc balloon	30 ml/cc balloon
12	MM51211210	MM51211230
14	MM51211410	MM51211430
16	MM51211610	MM51211630
18	MM51211810	MM51211830
20	MM51212010	MM51212030
22	MM51212210	MM51212230
24	MM51212410	MM51212430

**Material specifications:**

Contains latex? . . . . . Yes  
 Contains DEHP? . . . . . No  
 Catheter . . . . . Latex  
 Hard valve . . . . . Rubber  
 Sleeve . . . . . Rubber  
 Catheter surface . . . . . Silicone emulsion

**Certification:**

CE-marking? . . . . . Yes (CE 0123)  
 Product classification . . . . . Class IIB (EU)  
 Quality system certified to . . . . . ISO 13485  
 Global Medical Device Nomenclature (GMDN). . . . .34917

**Sterilisation:**

Sterile? . . . . . Yes  
 Sterilized by . . . . . Gamma radiation  
 Shelf life . . . . . 5 years

# Siliconized Foley catheter

Tiemann/Coudé, 2-way, soft valve

## Packaging:

Siliconized Foley catheter, Tiemann/Coudé, 2-way, soft valve

REF*	Shelf box		Shipping carton				MOQ, pcs.
	Qty., pcs	Dimensions, mm	Qty., pcs	Dimensions, mm	Gross Wt., kg	Volume, m <sup>3</sup>	
MM51211210	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM51211230	10	506 x 97 x 35	800	576 x 495 x 522	23	0.16	800
MM51211410	10	507 x 97 x 35	800	577 x 495 x 522	23	0.16	800
MM51211430	10	508 x 97 x 35	800	578 x 495 x 522	23	0.16	800
MM51211610	10	509 x 97 x 35	800	579 x 495 x 522	23	0.16	800
MM51211630	10	510 x 97 x 35	800	580 x 495 x 522	23	0.16	800
MM51211810	10	511 x 97 x 35	800	581 x 495 x 522	23	0.16	800
MM51211830	10	512 x 97 x 35	800	582 x 495 x 522	23	0.16	800
MM51212010	10	513 x 97 x 35	800	583 x 495 x 522	23	0.16	800
MM51212030	10	514 x 97 x 35	800	584 x 495 x 522	23	0.16	800
MM51212210	10	515 x 97 x 35	800	585 x 495 x 522	23	0.16	800
MM51212230	10	516 x 97 x 35	800	586 x 495 x 522	23	0.16	800
MM51212410	10	517 x 97 x 35	800	587 x 495 x 522	23	0.16	800
MM51212430	10	518 x 97 x 35	800	588 x 495 x 522	23	0.16	800

\* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

## Storage and handling:



*With reservations for misprints and alterations in product specifications. Printed version not controlled.*

Manufacturer: **Unomedical Sdn Bhd**  
(15290-P)  
Address: Bakar Arang Industrial Estate  
08000 Sungai Petani, Kedah, Malaysia

European Representative: **Unomedical a/s**  
Address: Kongevejen 2  
3460 Birkerød, Denmark

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.  
The ConvaTec logo is a trademark of ConvaTec Inc.

