

Inner cannula sets

for Unomedical tracheostomy tubes

Intended use:

A **tracheostomy tube** is inserted through a surgical opening (stoma) of the trachea to secure an airway during breathing, ventilation or weaning. The inner cannulae are inserted and locked in place after the obturator is removed. It acts as a removable liner for the more permanent, outer tube. The inner cannula can be withdrawn for brief periods for cleaning while the outer tube remains in place to maintain the patency of the airway. The product is intended for single use only.



List of product range:

Inner cannula sets, for usage with Unomedical TTIC

2 inner cannulas (non-fenestrated)	2 inner cannulas (1 fenestrated, 1 non-fenestrated)	Size	Tube, Inner diameter, mm
REF*	REF*		
MM72701060	MM72702060	6.0	5.0
MM72701070	MM72702070	7.0	5.5
MM72701075	MM72702075	7.5	6.0
MM72701080	MM72702080	8.0	6.5
MM72701085	MM72702085	8.5	7.0
MM72701090	MM72702090	9.0	7.5
MM72701100	MM72702100	10.0	8.0

* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Material specifications:

Contains latex? No
 Contains DEHP? No
 Inner cannula Polypropylene

Certification:

CE-marking? Yes (CE 0123)
 Product classification Class IIB (EU)
 Quality system certified to ISO 13485
 Global Medical Device Nomenclature (GMDN) 35404

Sterilisation:

Sterile? Yes
 Sterilized by Ethylene oxide
 Shelf life 5 years

Packaging:

Inner cannula sets, for usage with Unomedical TTIC

REF*	Shelf box		Shipping carton				SPU**, pcs.	MOQ, pcs.
	Qty., pcs.	Dimensions, mm	Qty., pcs.	Dimensions, mm	Gross Wt., kg	Volume, m ³		
MM72701XXX*	4	150 x 105 x 125	20	330 x 315 x 150	1.3	0.02	2	100
MM72702XXX*	4	150 x 105 x 125	20	330 x 315 x 150	1.3	0.02	2	100

* XXX = size, e.g 060 = 6.0. Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

** Smallest purchasing unit. 1 piece (set) contains 2 inner cannulas.

Storage and handling:



Inner cannula sets

for Unomedical tracheostomy tubes

With reservations for misprints and alterations in product specifications. Printed version not controlled.

Manufacturer: **Unomedical Sdn Bhd**
(15290-P)
Address: Bakar Arang Industrial Estate
08000 Sungai Petani, Kedah, Malaysia

European Representative: **Unomedical a/s**
Address: Kongevejen 2
3460 Birkerød, Denmark

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.
The ConvaTec logo is a trademark of ConvaTec Inc.

