


Siliconized Foley catheter

Standard, 2-way, hard valve

Intended use:

The siliconized Foley catheter is a urinary catheter with a silicone fluid treated surface for easy passage through the urethral meatus. The siliconized catheter is intended for short-term drainage of the bladder. The product is for single use only.



List of product range:

Siliconized Foley catheter, Standard, 2-way, hard valve

CH size	5-10 ml/cc balloon	30 ml/cc balloon
12	MM51121210	MM51121230
14	MM51121410	MM51121430
16	MM51121610	MM51121630
18	MM51121810	MM51121830
20	MM51122010	MM51122030
22	MM51122210	MM51122230
24	MM51122410	MM51122430
26	MM51122610	MM51122630
28	MM51122810	MM51122830
30	MM51123010	MM51123030

Material specifications:

Contains latex? Yes
 Contains DEHP? Yes
 Catheter Latex
 Hard valve Polypropylene
 Sleeve Polyvinyl chloride
 Catheter surface Silicone emulsion

Certification:

CE-marking? Yes (CE 0123)
 Product classification Class IIB (EU)
 Quality system certified to ISO 13485
 Global Medical Device Nomenclature (GMDN). 34917

Sterilisation:

Sterile? Yes
 Sterilized by Gamma radiation
 Shelf life 5 years

Siliconized Foley catheter

Standard, 2-way, hard valve

Packaging:

Siliconized Foley catheter, Standard, 2-way, hard valve

REF*	Shelf box		Shipping carton				MOQ, pcs.
	Qty., pcs	Dimensions, mm	Qty., pcs	Dimensions, mm	Gross Wt., kg	Volume, m ³	
MM51121210	10	505 x 95 x 34	800	575 x 495 x 522	23	0.16	800
MM51121230	10	505 x 95 x 34	800	575 x 495 x 522	23	0.16	800
MM51121410	10	505 x 95 x 34	800	575 x 495 x 522	23	0.16	800
MM51121430	10	505 x 95 x 34	800	575 x 495 x 522	23	0.16	800
MM51121610	10	505 x 95 x 34	800	575 x 495 x 522	23	0.16	800
MM51121630	10	505 x 95 x 34	800	575 x 495 x 522	23	0.16	800
MM51121810	10	505 x 95 x 34	800	575 x 495 x 522	23	0.16	800
MM51121830	10	505 x 95 x 34	800	575 x 495 x 522	23	0.16	800
MM51122010	10	505 x 95 x 34	800	575 x 495 x 522	23	0.16	800
MM51122030	10	505 x 95 x 34	800	575 x 495 x 522	23	0.16	800
MM51122210	10	505 x 95 x 34	800	575 x 495 x 522	23	0.16	800
MM51122230	10	505 x 95 x 34	800	575 x 495 x 522	23	0.16	800
MM51122410	10	505 x 95 x 34	800	575 x 495 x 522	23	0.16	800
MM51122430	10	505 x 95 x 34	800	575 x 495 x 522	23	0.16	800
MM51122610	10	505 x 95 x 34	800	575 x 495 x 522	23	0.16	800
MM51122630	10	505 x 95 x 34	800	575 x 495 x 522	23	0.16	800
MM51122810	10	505 x 95 x 34	800	575 x 495 x 522	23	0.16	800
MM51122830	10	505 x 95 x 34	800	575 x 495 x 522	23	0.16	800
MM51123010	10	505 x 95 x 34	800	575 x 495 x 522	23	0.16	800
MM51123030	10	505 x 95 x 34	800	575 x 495 x 522	23	0.16	800

* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Storage and handling:



With reservations for misprints and alterations in product specifications. Printed version not controlled.

Manufacturer: **Unomedical Sdn Bhd**
(15290-P)
Address: Bakar Arang Industrial Estate
08000 Sungai Petani, Kedah, Malaysia

European Representative: **Unomedical a/s**
Address: Kongevejen 2
3460 Birkerød, Denmark

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.
The ConvaTec logo is a trademark of ConvaTec Inc.

