

# Careline™ E4

## Sterile urine bag

**Intended use:**

Urine drainage bag with measuring scale for collection and monitoring of urine output. Fitted with bottom outlet.  
The product is intended for single use only.



**List of product range:**

Sterile urine bag, 2000 ml, 100 cm tubing, Needle-free KombiKon™ sample port and non-return valve (NRV)

	REF*	Cross valve	Lever tap
Cross valve - Hospital pack	38 126 1S-xx	Yes	No
Lever tap - Hospital pack	38 127 1S-xx	No	Yes
Cross valve - Community pack	38 132 1S-xx	Yes	No
Lever tap - Community pack	38 133 1S-xx	No	Yes

**Material specifications:**

**Contains latex?** . . . . . No  
**Contains DEHP?** . . . . . Yes  
**Tube inlet/outlet, foil** . . . . . Polyvinyl chloride  
**Connector with sample port.** . . . . . Polyvinyl chloride, Polypropylene, High density polyethylene, Polyethylene, Poly(styrene-butadiene-styrene)  
**Cap** . . . . . Low density polyethylene, High density polyethylene  
**Cross valve.** . . . . . Polypropylene, High density polyethylene  
**Lever tap** . . . . . Polyvinyl chloride, Acrylonitrile-butadiene-styrene  
**Films** . . . . . Polyvinyl chloride

**Certification:**

**CE-marking?** . . . . . Yes (CE 0086)  
**Product classification** . . . . . Class I (EU)  
**Quality system certified to** . . . . . ISO 13485  
**Global Medical Device Nomenclature (GMDN)** . . . . . 58921

**Sterilisation:**

**Sterile?** . . . . . Yes  
**Sterilized by** . . . . . Ethylene oxide  
**Shelf life** . . . . . 5 years

**Packaging:**

REF*	Inner carton		Shipping carton				MOQ, pcs.
	Qty., pcs.	Dimensions, mm	Qty., pcs.	Dimensions, mm	Gross Wt., kg	Volume, m <sup>3</sup>	
38 126 1S-xx	10	210 x 190 x 150	100	792 x 218 x 392	10.8	0.067	2,000
38 127 1S-xx	10	210 x 190 x 150	100	792 x 218 x 392	11.3	0.067	2,000
38 132 1S-xx	10	210 x 190 x 150	100	792 x 218 x 392	10.8	0.067	2,000
38 133 1S-xx	10	210 x 190 x 150	100	792 x 218 x 392	11.3	0.067	2,000

\* Here noted as xx, reference codes contain two additional digits (10-50) representing a language code.  
 Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

**Storage and handling:**



*With reservations for misprints and alterations in product specifications. Printed version not controlled.*

Manufacturer: **ConvaTec Limited**  
 First Avenue  
 Address: Deeside Industrial Park  
 Deeside, Flintshire, CH5 2NU, UK

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.  
 The ConvaTec logo is a trademark of ConvaTec Inc.