

Passive Wound Drainage



Wound drainage system

Intended use:

Passive wound drainage is a closed system designed for passive post-operative wound drainage. The passive wound drainage system is connected to a spiral drain which is placed directly into the surgical wound. Spiral drains, with recessed oval eyes in spiral formation, help reducing the risk of tissue contact and growth through opposing eyes. Passive wound drainage is primarily used after abdominal surgery in areas with delicate tissue. The product is intended for single use only.



Product range:

Passive wound drainage system with spiral drain	Drain bag capacity, ml	Tube length mm	Inner/outer diameter of drain, mm 					
			2.0/3.3	2.4/4.7	3.5/6.0	4.3/6.6	5.3/7.9	6.3/9.2
			CH 10 REF*	CH 14 REF*	CH 18 REF*	CH 20 REF*	CH 24 REF*	CH 28 REF*
	700	900	4022518x	4022618x	4021818x	4022918x	4023018x	4023118x

* Here noted as x, reference codes contain an additional digit (1-5) representing a language code. Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Material specifications:

Contains latex? No
 Contains DEHP? No
 Drainage bag Polyvinyl chloride
 Drain and tube Polyvinyl chloride

Certifications:

CE-marking? Yes (CE 0086)
 Product classification Class IIA (EU)
 Quality system certified to ISO 13485
 Global Medical Device Nomenclature (GMDN) 35824

Sterilization:

Sterile? Yes
 Sterilized by Ethylene oxide
 Shelf life 5 years

Packaging specifications:

REF	Shipping carton				MOQ, pcs
	Qty, pcs	Dimensions, mm (LxWxD)	Gross Wt. kg	Volume m ³	
4022518x	50	392 x 306 x 375	1.34	0.045	500
4022618x	50	392 x 306 x 375	1.34	0.045	500
4021818x	50	392 x 306 x 375	1.87	0.045	500
4022918x	50	392 x 306 x 375	4.90	0.045	500
4023018x	50	392 x 306 x 375	5.15	0.045	500
4023118x	50	392 x 306 x 375	5.50	0.045	500

Storage and handling:



ConvaTec Limited

First Avenue
 Deeside Industrial Park
 Deeside, Flintshire, CH5 2NU, UK

Distributed by



ConvaTec Limited
 First Avenue, Deeside Industrial Park
 Deeside, Flintshire, CH5 2NU, UK
www.convatec.com

For further information on these products, please contact your local ConvaTec sales representative.

With reservations for misprints and alterations in product specifications. Printed version not controlled.

All trade marks and logos are the property of a ConvaTec group company.

© 2018 ConvaTec Inc. AP-018718-MM, ver. FEB 2018