

POWERFUL PROTECTION

Against Skin Complications and Leaks.


Let's face it, for a person living with an ostomy, leaks can ruin everything - including their skin. And living with painful, red, irritated skin shouldn't be the norm.

It's time for a change.

It's time to **Go Moldable** with ConvaTec Moldable Technology™.

NOW Natura™ Accordion Flange with ConvaTec Moldable Technology™

Makes pouch application much easier and more comfortable - and delivers proven skin protection¹ for your patients.



SECURE:

Elastic-like seal “rebounds” to fit the size and shape of your patient’s stoma, providing an incredibly secure, snug fit. That means fewer leaks, and healthier skin.¹





SIMPLE:

Just roll, mold and hold the adhesive to accommodate any stoma shape and size. No cutting, no stretching, no guesswork. Flange lifts easily to provide full finger-width space – ample room for coupling.



COMFORTABLE:

Flexible skin barrier moves with the body, and can be put in “down” position to create a flat, low profile. The Natura™ Accordion Flange is compatible with our state-of-the-art comfort pouches.



The proven power of ConvaTec Moldable Technology™

- Over 95% of people who started with ConvaTec Moldable Technology™ kept their skin healthy.¹
- Over 86% of people who switched to ConvaTec Moldable Technology™ improved their skin condition.¹

Ordering Information:

Product Code	Color Code	Description	Stoma Size	Flange Size	Qty Per Box	HCPGS
421033	●	Moldable Stomahesive™ Wafer 57 mm	13-22 mm	57 mm	10	XXXXX
421034	●	Moldable Stomahesive™ Wafer 57 mm	22-33 mm	57 mm	10	XXXXX
421035	●	Moldable Stomahesive™ Wafer 70 mm	33-45 mm	70 mm	10	XXXXX
421039	●	Moldable Durahesive™ Wafer 57 mm	13-22 mm	57 mm	10	XXXXX
421040	●	Moldable Durahesive™ Wafer 57 mm	22-33 mm	57 mm	10	XXXXX
421041	●	Moldable Durahesive™ Wafer 70 mm	33-45 mm	70 mm	10	XXXXX

Find out More

1-800-422-8811 Monday – Friday 8:30 am – 7:00 pm

www.ConvaTec.com

References: 1. Szewczyk MT, Majewska GM, Cabral MV, Hölzel-Piontek K. Osmose Study: Multinational Evaluation of the Peristomal Condition in Ostomates Using Moldable Skin Barriers. Poster presented at ECET, Paris, France, June 2013.

™ indicates trademarks of ConvaTec Inc. Natura, Stomahesive and Durahesive are registered trademarks in the U.S.

© 2014 ConvaTec Inc. AP-014192-MM

