

Central-Gard™

Securement device

Intended use:

Central-Gard™ securement device is designed to efficiently secure subclavian and jugular lines. The product is intended for single use only.



List of product range:

	REF*	Size, cm
Central-Gard™, Small	667M	10.1 x 6.5
Central-Gard™, Large	668M	10.1 x 8.3

Material specifications:

Contains latex? No
Contains DEHP? No
Border and strips. Polyethylene foam
Adhesive. Acrylic adhesive
Film Polyurethane
Permeability Min. 1500g/m²/24h/37°C
Backing liner Polyethylene terephthalate
Paper Paper, Self-reinforcing polyethylene

Certification:

CE-marking? Yes (0086)
Product classification Class I (EU)
Quality system certified to ISO 13485
Global Medical Device Nomenclature (GMDN) 56631

Sterilisation:

Sterile? Yes
Sterilized by Gamma radiation
Shelf life 2 years

Packaging:

	REF*	Shelf box		Shipping carton				MOQ, pcs.
		Qty., pcs	Dimensions, mm	Qty., pcs	Dimensions, mm	Gross Wt., kg	Volume, m ³	
Central-Gard™, Small	667M	25	135 x 72 x 204	250	394 x 280 x 210	2.68	0.023	2,500
Central-Gard™, Large	668M	25	135 x 72 x 204	250	394 x 280 x 210	2.87	0.023	2,500

* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Storage and handling:



With reservations for misprints and alterations in product specifications. Printed version not controlled.

Manufacturer: **ConvaTec Limited**
 First Avenue
 Address: Deeside Industrial Park
 Deeside, Flintshire, CH5 2NU, UK

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.
 The ConvaTec logo is a trademark of ConvaTec Inc.