

Endotracheal Tube

Standard & Preformed

Intended use:

The Endotracheal Tube (ETT) is inserted through the mouth or nose into the trachea and connected to a ventilation device to secure an airway. The product is intended for single use only.



List of product range:

Standard Endotracheal Tubes with cuff

Magill shaped ETT with HVLP* cuff		Magill shaped ETT with EVLP* cuff						
With Murphy eye	Without Murphy eye	With Murphy eye	Without Murphy eye					
REF	REF	REF	REF	Size	I.D. mm	O.D. mm	Length** mm	Connector color
MM61110030	MM61120030			3.0	3.0	4.1	185	Red
MM61110035	MM61120035			3.5	3.5	4.8	185	Red
MM61110040	MM61120040			4.0	4.0	5.4	205	Green
MM61110045	MM61120045			4.5	4.5	6.1	225	Green
MM61110050	MM61120050			5.0	5.0	6.7	245	Blue
MM61110055	MM61120055			5.5	5.5	7.3	275	Blue
MM61110060	MM61120060	MM61410060	MM61420060	6.0	6.0	8.0	285	Blue
MM61110065	MM61120065	MM61410065	MM61420065	6.5	6.5	8.7	295	Black
MM61110070	MM61120070	MM61410070	MM61420070	7.0	7.0	9.3	305	Black
MM61110075	MM61120075	MM61410075	MM61420075	7.5	7.5	10.0	315	White
MM61110080	MM61120080	MM61410080	MM61420080	8.0	8.0	10.7	325	White
MM61110085	MM61120085	MM61410085	MM61420085	8.5	8.5	11.3	325	Green
MM61110090	MM61120090	MM61410090	MM61420090	9.0	9.0	12.0	325	Green
MM61110095	MM61120095			9.5	9.5	12.7	325	Orange
MM61110100	MM61120100			10.0	10.0	13.3	325	Orange

* HVLP: High Volume Low Pressure, EVLP: Extended Volume Low Pressure

** Length is measured from tip to angle

Endotracheal Tube

Standard & Preformed

List of product range:

Standard Endotracheal Tubes without cuff

Magill shaped ETT without cuff						
With Murphy eye	Without Murphy eye					
REF	REF	Size	I.D. mm	O.D. mm	Length**, mm	Connector color
MM61130020	MM61140020	2.0	2.0	2.9	150	Clear
MM61130025	MM61140025	2.5	2.5	3.4	175	Clear
MM61130030	MM61140030	3.0	3.0	4.1	185	Red
MM61130035	MM61140035	3.5	3.5	4.8	185	Red
MM61130040	MM61140040	4.0	4.0	5.4	205	Green
MM61130045	MM61140045	4.5	4.5	6.1	225	Green
MM61130050	MM61140050	5.0	5.0	6.7	245	Blue
MM61130055	MM61140055	5.5	5.5	7.3	275	Blue
MM61130060	MM61140060	6.0	6.0	8.0	285	Blue
MM61130065	MM61140065	6.5	6.5	8.7	295	Black
MM61130070	MM61140070	7.0	7.0	9.3	305	Black
MM61130075	MM61140075	7.5	7.5	10.0	315	White
MM61130080	MM61140080	8.0	8.0	10.7	325	White
MM61130085	MM61140085	8.5	8.5	11.3	325	Green
MM61130090	MM61140090	9.0	9.0	12.0	325	Green
MM61130095	MM61140095	9.5	9.5	12.0	325	Orange
MM61130100	MM61140100	10.0	10.0	13.3	325	Orange

List of product range:

Oral preformed

With HVLP* cuff with Murphy eye	Without cuff with Murphy eye					
REF	REF	Size	I.D. mm	O.D. mm	Length**, mm	Connector color
	MM61233030	3.0	3.0	4.1	105	Red
	MM61233035	3.5	3.5	4.8	120	Red
MM61213040	MM61233040	4.0	4.0	5.4	130	Green
MM61213045	MM61233045	4.5	4.5	6.1	150	Green
MM61213050	MM61233050	5.0	5.0	6.7	155	Blue
MM61213055	MM61233055	5.5	5.5	7.3	165	Blue
MM61213060	MM61233060	6.0	6.0	8.0	180	Blue
MM61213065	MM61233065	6.5	6.5	8.7	190	Black
MM61213070	MM61233070	7.0	7.0	9.3	200	Black
MM61213075		7.5	7.5	10.0	210	White
MM61213080		8.0	8.0	10.7	220	White
MM61213085		8.5	8.5	11.3	230	Green
MM61213090		9.0	9.0	11.3	240	Green

* HVLP: High Volume Low Pressure

** Length is measured from tip to angle

Endotracheal Tube

Standard & Preformed

List of product range:

Nasal preformed

With HVLP* cuff with Murphy eye REF	With HVLP* cuff without Murphy eye REF	Without cuff with Murphy eye REF	Size	I.D., mm	O.D., mm	Without cuff Length**, mm	With cuff Length**, mm	Connector color
		MM61333030	3.0	3.0	4.1	131		Red
		MM61333035	3.5	3.5	4.8	160		Red
		MM61333040	4.0	4.0	5.4	175		Green
		MM61333045	4.5	4.5	6.1	202		Green
		MM61333050	5.0	5.0	6.7	212		Blue
		MM61333055	5.5	5.5	7.3	215		Blue
MM61313060	MM61323060	MM61333060	6.0	6.0	8.0	227	295	Blue
MM61313065	MM61323065	MM61333065	6.5	6.5	8.7	234	305	Black
MM61313070	MM61323070	MM61333070	7.0	7.0	9.3	237	313	Black
MM61313075	MM61323075		7.5	7.5	10.0		328	White
MM61313080	MM61323080		8.0	8.0	10.7		335	White

List of product range:

Oxford tube

With HVLP* cuff REF	Size	I.D. mm	O.D. mm	Length**, mm	Connector color
MM61222060	6.0	6.0	8.0	195	Blue
MM61222070	7.0	7.0	9.3	215	Black
MM61222080	8.0	8.0	10.7	225	White
MM61222090	9.0	9.0	12.0	235	Green

* HVLP: High Volume Low Pressure

** Length is measured from tip to angle

Material specifications:

Contains latex? No
 Contains DEHP? Yes
 Tube Polyvinyl chloride
 Cuff Polyvinyl chloride
 Connector Polypropylene
 Pilot balloon Polyvinyl chloride
 Inflation line Polyvinyl chloride
 Valve Polyvinyl chloride, Nitrile, Stainless steel

Certification:

CE-marking? Yes (CE 0123)
 Product classification Class IIA (EU)
 Quality system certified to ISO 13485
 Global Medical Device Nomenclature (GMDN) 46967

Sterilisation:

Sterile? Yes
 Sterilized by Ethylene oxide
 Shelf life 5 years

Endotracheal Tube

Standard & Preformed

Packaging:

ETT standard models

REF*	Size	Shelf box		Shipping carton				MOQ, pcs.
		Qty., pcs	Dimensions, mm	Qty., pcs	Dimensions, mm	Gross Wt., kg	Volume, m ³	
MM61110XXX	3.0 to 10.0	10	348 x 150 x 98	100	540 x 380 x 340	6.2	0.07	300
MM61120XXX	3.0 to 10.0	10	348 x 150 x 98	100	540 x 380 x 340	6.2	0.07	300
MM61130XXX	2.0 to 10.0	10	348 x 150 x 98	100	540 x 380 x 340	6.2	0.07	300
MM61140XXX	2.0 to 10.0	10	348 x 150 x 98	100	540 x 380 x 340	6.2	0.07	300
MM61410XXX	6.0 to 9.0	10	348 x 150 x 98	100	540 x 380 x 340	6.2	0.07	300
MM61420XXX	6.0 to 9.0	10	348 x 150 x 98	100	540 x 380 x 340	6.2	0.07	300
MM63130XXX	2.0 to 6.0	10	348 x 150 x 98	100	540 x 380 x 340	6.2	0.07	300
MM63140XXX	2.0 to 6.0	10	348 x 150 x 98	100	540 x 380 x 340	6.2	0.07	300

Packaging:

Preformed ETT models

REF*	Size	Shelf box		Shipping carton				MOQ, pcs.
		Qty., pcs	Dimensions, mm	Qty., pcs	Dimensions, mm	Gross Wt., kg	Volume, m ³	
MM61213XXX	4.0 to 9.0	10	348 x 150 x 98	100	540 x 380 x 340	6.2	0.07	300
MM61233XXX	3.0 to 7.0	10	348 x 150 x 98	100	540 x 380 x 340	6.2	0.07	300
MM61313XXX	6.0 to 8.0	10	388 x 202 x 98	100	525 x 425 x 415	6.8	0.07	300
MM61323XXX	6.0 to 8.0	10	388 x 202 x 98	100	525 x 425 x 415	6.8	0.07	300
MM61333XXX	3.0 to 7.0	10	348 x 150 x 98	100	540 x 380 x 340	6.2	0.07	300
MM61222XXX	6.0, 7.0, 8.0, 9.0	10	348 x 150 x 98	100	540 x 380 x 340	6.2	0.07	300

*XXX= size, e.g. 045 = 4.5 mm I.D.

* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Storage and handling:



With reservations for misprints and alterations in product specifications. Printed version not controlled.

Manufacturer: **Unomedical Sdn Bhd**
(15290-P)
Address: Bakar Arang Industrial Estate
08000 Sungai Petani, Kedah, Malaysia

European Representative: **Unomedical a/s**
Address: Kongevejen 2
3460 Birkerød, Denmark

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.
The ConvaTec logo is a trademark of ConvaTec Inc.

