

ClearVue™ Anesthesia Mask















with pre-inflated/adjustable cushion

Intended use:

Anesthesia masks (face masks) are intended for the delivery of anesthetic gases, air and/or oxygen to patients. The products are for single use only.



Product range:

	REF*	Item	Size		Hook ring (removable)	Air valve
			Mask	Connector		
	93-580UM	ClearVue™ Anesthesia Mask with pre-inflated cushion	Size 0 Premature	15mm (Male)	-	-
	93-581UM	ClearVue™ Anesthesia Mask with pre-inflated cushion	Size 1 Neonate	15mm (Male)	-	-
	93-582UM	ClearVue™ Anesthesia Mask with pre-inflated cushion	Size 2 Infant	15mm (Male)	-	-
	93-583UM	ClearVue™ Anesthesia Mask with adjustable cushion	Size 3 Pediatric	22mm (Female)	White	Yes
	93-584UM	ClearVue™ Anesthesia Mask with adjustable cushion	Size 4 Pediatric/Adult (Small)	22mm (Female)	Yellow	Yes
	93-585UM	ClearVue™ Anesthesia Mask with adjustable cushion	Size 5 Adult (Medium)	22mm (Female)	Green	Yes
	93-586UM	ClearVue™ Anesthesia Mask with adjustable cushion	Size 6 Adult (Large)	22mm (Female)	Orange	Yes
	93-580SM (Strawberry Scented)	ClearVue™ Anesthesia Mask with pre-inflated cushion, Scented	Size 0 Premature	15mm (Male)	-	-
	93-581SM (Strawberry Scented)	ClearVue™ Anesthesia Mask with pre-inflated cushion, Scented	Size 1 Neonate	15mm (Male)	-	-
	93-582SM (Strawberry Scented)	ClearVue™ Anesthesia Mask with pre-inflated cushion, Scented	Size 2 Infant	15mm (Male)	-	-
	93-583SM (Mint Scented)	ClearVue™ Anesthesia Mask with adjustable cushion, Scented	Size 3 Pediatric	22mm (Female)	White	Yes
	93-584SM (Mint Scented)	ClearVue™ Anesthesia Mask with adjustable cushion, Scented	Size 4 Pediatric/Adult (Small)	22mm (Female)	Yellow	Yes
	93-585SM (Mint Scented)	ClearVue™ Anesthesia Mask with adjustable cushion, Scented	Size 5 Adult (Medium)	22mm (Female)	Green	Yes
	93-586SM (Mint Scented)	ClearVue™ Anesthesia Mask with adjustable cushion, Scented	Size 6 Adult (Large)	22mm (Female)	Orange	Yes

*Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

ClearVue™ Anesthesia Mask

with pre-inflated/adjustable cushion



Material specifications:

Contains latex?	No
Contains phthalates?	No
Contains bisphenols?	No
Anesthesia mask	Polycarbonate
Cushion	Polyvinyl chloride
Connector	Polycarbonate
Hook ring	Polycarbonate
Air valve	Stainless steel, Silicone, Polyvinyl chloride Acrylonitrile-butadiene-styrene

Certifications:

CE-marking?	Yes (CE 1639)
Product classification	Class IIA (EU)
Quality system certified to	ISO 13485

Sterilization:

Sterile?	No
Shelf life	3 years

Packaging specifications:

REF	Inner package		Outer package/Shipper				Pallet
	Qty, pcs	Dimensions, mm (LxW)	Qty, pcs	Dimensions, mm (LxWxD)	Gross Wt. kg	Volume m ³	Qty, pcs
93-580UM	1	170 x 180	20	280 x 200 x 175	0.60	0.010	1,920
93-581UM	1	170 x 180	20	280 x 200 x 175	0.65	0.010	1,920
93-582UM	1	170 x 180	20	280 x 200 x 175	0.75	0.010	1,920
93-583UM	1	190 x 180	50	400 x 300 x 240	2.40	0.029	1,600
93-584UM	1	210 x 180	50	400 x 380 x 260	2.90	0.040	1,200
93-585UM	1	180 x 235	50	400 x 400 x 350	3.80	0.056	900
93-586UM	1	200 x 255	50	600 x 400 x 260	4.80	0.062	800
93-580SM	1	170 x 180	20	280 x 200 x 175	0.60	0.010	1,920
93-581SM	1	170 x 180	20	280 x 200 x 175	0.65	0.010	1,920
93-582SM	1	170 x 180	20	280 x 200 x 175	0.75	0.010	1,920
93-583SM	1	190 x 180	50	400 x 300 x 240	2.40	0.029	1,600
93-584SM	1	210 x 180	50	400 x 380 x 260	2.90	0.040	1,200
93-585SM	1	180 x 235	50	400 x 400 x 350	3.80	0.056	900
93-586SM	1	200 x 255	50	600 x 400 x 260	4.80	0.062	800

EAN codes:

REF	Inner package	Outer package
93-580UM	(01)06945099304278	(01)46945099304276
93-581UM	(01)06945099304285	(01)46945099304283
93-582UM	(01)06945099304292	(01)46945099304290
93-583UM	(01)06945099304308	(01)46945099304306
93-584UM	(01)06945099304315	(01)46945099304313
93-585UM	(01)06945099304322	(01)46945099304320
93-586UM	(01)06945099304339	(01)46945099304337
93-580SM	(01)06945099304346	(01)46945099304344
93-581SM	(01)06945099304353	(01)46945099304351
93-582SM	(01)06945099304360	(01)46945099304368
93-583SM	(01)06945099304377	(01)46945099304375
93-584SM	(01)06945099304384	(01)46945099304382
93-585SM	(01)06945099304391	(01)46945099304399
93-586SM	(01)06945099304407	(01)46945099304405



ClearVue™ Anesthesia Mask





with pre-inflated/adjustable cushion



Packaging material:

Inner package	Poly bag - Low density polyethylene / Insert - 80g paper
Outer package	Carton - Corrugated board / Sticker - Coated paper

Storage and handling:

	FRAGILE - Contents of the transport package are fragile and shall therefore be handled with care
	KEEP DRY
	THIS WAY UP - Indicates correct upright position of the transport package
	DO NOT USE HAND HOOKS - Hooks are prohibited for handling the transport package



Ningbo Shengyurui Medical Appliances Co., Ltd.

No. 138 Binhaisi Rd.
Hangzhou Bay New Zone, 315336, Ningbo City
Zhejiang Province, P.R. China

Made in China



Shanghai Intl. Holding Corp. GmbH

Eiffestrasse 80, 20537 Hamburg
Germany

Distributed by



ConvaTec Limited
First Avenue, Deeside Industrial Park
Deeside, Flintshire, CH5 2NU, UK

www.convatec.com

For further information on these products, please contact your local ConvaTec sales representative.

With reservations for misprints and alterations in product specifications. Printed version not controlled.
All trade marks and logos are the property of a ConvaTec group company.

© 2020 ConvaTec Inc. AP-021052-MM, ver. AUG 2020