

SimpaVac™

Wound drainage system

Intended use:

SimpaVac™ is a low vacuum wound drainage system designed for post-operative wound drainage procedure. The product is intended for single use only.



List of product range:

REF*	REF*
SimpaVac™, 200 ml	4007118x

Material specifications:

Contains latex? No
Contains DEHP? No
O-ring Silicone
Container Polyethylene
Straps Polyethylene
Plug Polycarbonate

Certification:

CE-marking? Yes (CE 0086)
Product classification Class 1 sterile (EU)
Quality system certified to ISO 13485
Global Medical Device Nomenclature (GMDN) 35824

Sterilisation:

Sterile? Yes
Sterilized by Ethylene oxide
Shelf life 5 years

Packaging:

REF*	Shipping carton				MOQ, pcs.
	Qty., pcs.	Dimensions, mm	Gross Wt., kg	Volume, m ³	
4007118x	50	391 x 304 x 455	3.2	0.056	100

* Here noted as x, reference codes contain an additional digit (1-5) representing a language code.
 Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Storage and handling:



With reservations for misprints and alterations in product specifications. Printed version not controlled.

Manufacturer: **ConvaTec Limited**
 First Avenue
 Address: Deeside Industrial Park
 Deeside, Flintshire, CH5 2NU, UK

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.
 The ConvaTec logo is a trademark of ConvaTec Inc.