


HandyVac™


Wound drainage system

Intended use:

Low vacuum wound drainage system for post-operative wound drainage. The product is intended for single use only.



Product range:

	Drain				Tube		Bellows	Collection Bag				Inner/outer diameter of drain, mm 					
	Spiral drain	Drain length, mm	Perforation length, mm	Trocar	Length, mm	Connector to drain removable	Capacity, ml	Capacity, ml	Bag inlet with closable cap	Bottom outlet	1.6/2.7	2.0/3.3	2.2/4.0	2.4/4.7	3.0/5.3	3.5/6.0	
											CH 08 REF*	CH 10 REF*	CH 12 REF*	CH 14 REF*	CH 16 REF*	CH 18 REF*	
HandyVac™ 700 ml bag with bottom outlet	Yes	650	150	-	800	Yes	125	700	-	Yes	4009918x	4010118x	4010318x	4010518x	4017118x	4010718x	
	Yes	650	150	Yes	800	Yes	125	700	-	Yes	4012118x	4012318x	4016018x	4012518x	4018418x	4012718x	
HandyVac™ 700 ml bag without bottom outlet	Yes	650	150	Yes	800	Yes	125	700	Yes	-	4022118x	4022218x		4022318x		4022418x	
	Yes	650	150	Yes	800	Yes	125	700	-	-			4022718x		4022818x		
	Yes	650	150	Yes	800	-	125	700	Yes	-	4023218x	4023318x	4023418x	4023518x	4023618x	4023718x	
HandyVac™ 700 ml bag with bottom outlet, no drain	-	-	-	-	800	-	125	700	-	Yes						4009718x	
HandyVac™ 700 ml bag without bottom outlet, no drain	-	-	-	-	800	-	125	700	-	-						4026618x	
Extra drainage bag 700 ml with bottom outlet	-	-	-	-	-	-	-	700	-	Yes						40067069	
Extra drainage bag 700 ml with bottom outlet	-	-	-	-	900	-	-	700	-	Yes						40069069	
Extra drainage bag 700 ml without bottom outlet	-	-	-	-	-	-	-	700	Yes	-						4020018x	

* Here noted as x, reference codes contain an additional digit (1-5) representing a language code. Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Material specifications:

- Contains latex? No
- Contains DEHP? No
- Spiral drain Polyvinyl chloride
- Trocar Stainless steel
- Tube clamp High density polyethylene
- HandyVac™ bellows Ethylene butyl acrylate
- Strap Low density polyethylene
- Strap for plug Ethylene butyl acrylate
- Plug Polycarbonate
- Drainage bag Polyvinyl chloride, Polycarbonate, High density polyethylene, Expanded polytetrafluoroethylene, Polypropylene
- Bushings Polyvinyl chloride
- Connectors Polycarbonate
- Swivel for connector Copolyester
- O-ring Silicone
- Foil spreader Polyvinyl chloride

Certifications:

- CE-marking? Yes (CE 0086)
- Product classification Class I (EU) (without drain)
Class IIA (EU) (with drain)
- Quality system certified to ISO 13485
- Global Medical Device Nomenclature (GMDN) 35824

Sterilization:

- Sterile? Yes
- Sterilized by Ethylene oxide
- Shelf life 5 years

Packaging specifications:

REF	Shelf box		Shipping carton				MOQ, pcs
	Qty, pcs	Dimensions, mm (LxWxD)	Qty, pcs	Dimensions, mm (LxWxD)	Gross Wt. kg	Volume m ³	
4009918x	*	*	20	392 x 306 x 375	3.20	0.045	100
4010118x	*	*	20	392 x 306 x 375	3.20	0.045	100
4010318x	*	*	20	392 x 306 x 375	3.20	0.045	100
4010518x	*	*	20	392 x 306 x 375	3.20	0.045	100
4017118x	*	*	20	392 x 306 x 375	3.50	0.045	60
4010718x	*	*	20	392 x 306 x 375	3.50	0.045	100
4012118x	*	*	20	392 x 306 x 375	3.50	0.045	100
4012318x	*	*	20	392 x 306 x 375	3.60	0.045	100
4016018x	*	*	20	392 x 306 x 375	4.00	0.045	100
4012518x	*	*	20	392 x 306 x 375	4.00	0.045	100
4018418x	*	*	20	392 x 306 x 375	4.02	0.045	20
4012718x	*	*	20	392 x 306 x 375	4.15	0.045	100
4022118x	*	*	20	392 x 306 x 375	3.50	0.045	100
4022218x	*	*	20	392 x 306 x 375	3.60	0.045	100
4022318x	*	*	20	392 x 306 x 375	4.00	0.045	100
4022418x	*	*	20	392 x 306 x 375	4.15	0.045	100
4022718x	*	*	20	392 x 306 x 375	3.60	0.045	100
4022818x	*	*	20	392 x 306 x 375	4.02	0.045	100
4023218x	*	*	20	392 x 306 x 375	3.35	0.045	100
4023318x	*	*	20	392 x 306 x 375	3.40	0.045	100
4023418x	*	*	20	392 x 306 x 375	3.40	0.045	100
4023518x	*	*	20	392 x 306 x 375	4.00	0.045	100
4023618x	*	*	20	392 x 306 x 375	4.00	0.045	100
4023718x	*	*	20	392 x 306 x 375	4.15	0.045	100
4009718x	*	*	20	392 x 306 x 375	3.00	0.045	100
4026618x	*	*	20	392 x 306 x 375	3.00	0.045	100
40067069	100	580 x 380 x 176	200	591 x 391 x 385	8.10	0.089	200
40069069	100	580 x 380 x 176	200	591 x 391 x 385	14.00	0.089	200
4020018x	50	294 x 125 x 360	150	392 x 306 x 375	13.30	0.045	150

* Single packed, No shelf box

Storage and handling:




ConvaTec Limited

First Avenue
Deeside Industrial Park
Deeside, Flintshire, CH5 2NU, UK

Distributed by



ConvaTec Limited
First Avenue, Deeside Industrial Park
Deeside, Flintshire, CH5 2NU, UK

www.convatec.com

For further information on these products, please contact your local ConvaTec sales representative.

With reservations for misprints and alterations in product specifications. Printed version not controlled.

All trade marks and logos are the property of a ConvaTec group company.

© 2018 ConvaTec Inc. AP-018717-MM, ver. FEB 2018