

Careline™ Combi Bag

Urine drainage bag

Intended use:

Urine drainage bag with scale for collection and monitoring of urine output. The device is not intended to have a measurement function. The scale is approximate and for indication purposes only.



List of product range:

REF*	Volume, ml	Tube length, cm	Cross valve	Cut to fit
32 450 0-xx	1,500	90	Yes	Yes

Material specifications:

Contains latex? No
Contains DEHP? Yes
Tube Polyurethane with X-ray line
Connector Polyvinyl chloride
Cap for connector Low density polyethylene,
 High density polyethylene
Adjustable connector Polypropylene,
 High density polyethylene
Cap for adjustable connector Low density polyethylene,
 High density polyethylene
Cross valve Polypropylene, High density polyethylene

Certification:

CE-marking? Yes
Product classification Class I (EU)
Quality system certified to ISO 13485
Global Medical Device Nomenclature (GMDN) 58922

Sterilisation:

Sterile? No
Shelf life 5 years

Packaging:

REF*	Shipping carton				MOQ, pcs.
	Qty, pcs.	Dimensions, mm	Gross Wt., kg	Volume, m ³	
32 450 0-xx	30	375 x 280 x 81	2.4	0.009	3,840

* Here noted as xx, reference codes contain two additional digits (10-50) representing a language code.
 Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Storage and handling:



With reservations for misprints and alterations in product specifications. Printed version not controlled.

Manufacturer: **ConvaTec Limited**
 First Avenue
 Address: Deeside Industrial Park
 Deeside, Flintshire, CH5 2NU, UK

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.
 The ConvaTec logo is a trademark of ConvaTec Inc.