


# Sump Yankauer

## Suction handles and sets

**Intended use:**

The Sump Yankauer is a suction handle intended for the removal of blood and fluids during surgery. The Sump Yankauer is designed with a diffuser tip for delicate tissue. The product is intended for single use only.



**List of product range:**

	REF*	Connecting tube Shore A 74		Inner/outer diameter of suction, mm	Suction handle, length in mm	Vacuum control	Number of relieving eyes	CH/FG (handle)	
		Length in mm	CH/FG						
Sump Yankauer, diffusor tip, single bend	3301518x	-	-	4.6/6.5	210	No	0	20	Medium
	3401018x	2,100	24	4.6/6.5	210	No	0	20	Medium
	3402118x	2,100	24	4.6/6.5	210	Yes	0	20	Medium

**Material specifications:**

**Contains latex?** . . . . . No  
**Contains DEHP** . . . . . Yes  
**Handle & tip** . . . . . Polyvinyl chloride  
**Tubing** . . . . . Polyvinyl chloride  
**Connector** . . . . . Polyvinyl chloride

**Certification:**

**CE-marking?** . . . . . Yes (CE 0086)  
**Product classification** . . . . . Class IIA (EU)  
**Quality system certified to** . . . . . ISO 13485  
**Global Medical Device**  
**Nomenclature (GMDN)** . . . . .  
 (without connecting tube) 35917, (with connecting tube) 36778

**Sterilisation:**

**Sterile?** . . . . . Yes  
**Sterilized by** . . . . . Ethylene oxide  
**Shelf life** . . . . . 5 years

**Packaging:**

REF*	Shelf box		Shipping carton				MOQ, pcs.
	Qty., pcs	Dimensions, mm	Qty., pcs	Dimensions, mm	Gross Wt., kg	Volume, m <sup>3</sup>	
3301518x	50	294 x 125 x 360	300	391 x 304 x 743	5.90	0.092	600
3401018x	*	*	50	392 x 306 x 340	6.00	0.043	600
3402118x	*	*	50	392 x 306 x 340	6.00	0.043	600

\* = Single packed - no shelf box

\* Here noted as x, reference codes contain an additional digit (1-5) representing a language code.  
 Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

**Storage and handling:**



*With reservations for misprints and alterations in product specifications. Printed version not controlled.*

Manufacturer: **ConvaTec Limited**  
 First Avenue  
 Address: Deeside Industrial Park  
 Deeside, Flintshire, CH5 2NU, UK

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.  
 The ConvaTec logo is a trademark of ConvaTec Inc.