

Careline™ All Silicone Foley catheter

Standard, 2-way, hard valve, with prefilled syringe

Intended use:

Urethral Foley catheter for drainage of urine from the bladder. Must be connected to a drainage device.
The product is intended for single use only.



List of product range:

Careline™ All Silicone Foley catheter, Standard, 2 way, hard valve, with prefilled syringe

CH size	5-10 ml/cc balloon
12	MM41151205P
14	MM41151405P
16	MM41151605P
18	MM41151805P
20	MM41152005P
22	MM41152205P
24	MM41152405P

Material specifications:

Contains latex? No
Contains DEHP? Yes
Catheter tube and balloon Silicone
Valve Polypropylene
Sleeve Polyvinyl chloride
Tip Silicone, Barium sulphate
Syringe Polypropylene, Latex free elastomer
Pre-filled content. Sterile water

Certification:

CE-marking? Yes (CE 0123)
Product classification Class IIB (EU)
Quality system certified to ISO 13485
Global Medical Device Nomenclature (GMDN) 34917

Sterilisation:

Sterile? Yes
Sterilized by Gamma radiation
Shelf life 2 years

Packaging:

Careline™ All Silicone Foley catheter, Standard, 2 way, hard valve, with prefilled syringe

REF*	Shelf box		Shipping carton				MOQ. pcs.
	Qty. pcs	Dimensions, mm	Qty. pcs	Dimensions, mm	Gross Wt., kg	Volume, m ³	
MM41151205P	5	500 x 92 x 52	250	580 x 500 x 527	9.50	0.15	500
MM41151405P	5	500 x 92 x 52	250	580 x 500 x 527	10.00	0.15	500
MM41151605P	5	500 x 92 x 52	250	580 x 500 x 527	13.00	0.15	500
MM41151805P	5	500 x 92 x 52	250	580 x 500 x 527	14.00	0.15	500
MM41152005P	5	500 x 92 x 52	250	580 x 500 x 527	15.00	0.15	500
MM41152205P	5	500 x 92 x 52	250	580 x 500 x 527	16.00	0.15	500
MM41152405P	5	500 x 92 x 52	250	580 x 500 x 527	17.00	0.15	500

* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Storage and handling:



Careline™ All Silicone Foley catheter

Standard, 2-way, hard valve, with prefilled syringe

With reservations for misprints and alterations in product specifications. Printed version not controlled.

Manufacturer: **Unomedical Sdn Bhd**
(15290-P)
Address: Bakar Arang Industrial Estate
08000 Sungai Petani, Kedah, Malaysia

European Representative: **Unomedical a/s**
Address: Kongevejen 2
3460 Birkerød, Denmark

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.
The ConvaTec logo is a trademark of ConvaTec Inc.

