


Pediatric Foley catheters

Coated Latex

Intended use:

The pediatric foley catheter is a urinary catheter intended for short as well as long term drainage of the bladder. For routine and post-operative drainage. The size is specially designed to fit pediatric patients. The product is intended for single use only.


List of product range:

Pediatric Foley catheter, Latex

Siliconized Foley catheter, Pediatric, Standard, 2-way, hard valve

Balloon size	CH 06	CH 08	CH 10
3-5 ml/cc	MM53120605	MM53120805	MM53121005

Silicone elastomer coated Foley catheter, Pediatric, Standard, 2-way, hard valve

Balloon size	CH 06	CH 08	CH 10
3-5 ml/cc	MM53160605	MM53160805	MM53161005

Lubricious coated Foley catheter, Pediatric, Standard, 2-way, hard valve

Balloon size	CH 06	CH 08	CH 10
3-5 ml/cc	MM53180605	MM53180805	MM53181005

Material specifications:
Siliconized Foley catheter

Contains latex? Yes
 Contains DEHP? Yes
 Catheter tube and balloon Latex
 Hard valve Polypropylene
 Sleeve Polyvinyl chloride
 Catheter surface Silicone emulsion

Silicone Elastomer coated Foley catheter

Contains latex? Yes
 Contains DEHP? Yes
 Catheter tube and balloon Latex
 Hard valve Polypropylene
 Sleeve Polyvinyl chloride
 Catheter surface Silicone elastomer coating

Lubricious coated Foley catheter

Contains latex? Yes
 Contains DEHP? Yes
 Catheter tube and balloon Latex
 Hard valve Polypropylene
 Sleeve Polyvinyl chloride
 Catheter surface Hydrophilic coating

Certification:

CE-marking? Yes (CE 0123)
 Product classification Class IIB (EU)
 Quality system certified to ISO 13485
 Global Medical Device Nomenclature (GMDN) 34917

Sterilisation:

Sterile? Yes
 Sterilized by Gamma radiation
 Shelf life 5 years

Pediatric Foley catheters

Coated Latex

Packaging:

REF*	Shelf box		Shipping carton				MOQ, pcs.
	Qty., pcs	Dimensions, mm	Qty., pcs	Dimensions, mm	Gross Wt., kg	Volume, m ³	
MM53120605	10	505 x 97 x 35	800	575 x 495 x 522	23,00	0.16	800
MM53120805	10	505 x 97 x 35	800	575 x 495 x 522	23,00	0.16	800
MM53121005	10	505 x 97 x 35	800	575 x 495 x 522	23,00	0.16	800
MM53160605	10	505 x 97 x 35	800	575 x 495 x 522	23,00	0.16	800
MM53160805	10	505 x 97 x 35	800	575 x 495 x 522	23,00	0.16	800
MM53161005	10	505 x 97 x 35	800	575 x 495 x 522	23,00	0.16	800
MM53180605	10	505 x 97 x 35	800	575 x 495 x 522	23,00	0.16	800
MM53180805	10	505 x 97 x 35	800	575 x 495 x 522	23,00	0.16	800
MM53181005	10	505 x 97 x 35	800	575 x 495 x 522	23,00	0.16	800

* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Storage and handling:



With reservations for misprints and alterations in product specifications. Printed version not controlled.

Manufacturer: **Unomedical Sdn Bhd**
(15290-P)
Address: Bakar Arang Industrial Estate
08000 Sungai Petani, Kedah, Malaysia

European Representative: **Unomedical a/s**
Address: Kongevejen 2
3460 Birkerød, Denmark

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.
The ConvaTec logo is a trademark of ConvaTec Inc.

