


Bacterial/Viral Filters and Heat & Moisture Exchangers





Intended use:

Bacterial/Viral Filters and Heat & Moisture Exchangers provide filtration and moisture return required during anesthesia and critical care. The product is intended for single use only.



List of product range:

Connectors are standard
one end: ID 15 mm, OD 22 mm; other end: ID 22 mm, OD 25.5 mm

	REF*	Description
	761MM	Heat & moisture exchanger
	762MM	Heat & moisture exchanger with Luer lock sample port
	764MM	Heat & moisture exchanger with Bacterial/Viral filter and Luer lock sample port
	864MM	Bacterial/Viral Filter

* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Material specifications:

Contains latex? No
Contains DEHP? No
Heat & moisture exchanger Polyurethane
Bacterial/Viral filter K-Resin
Connectors Polypropylene
Luer lock sample port K-Resin

Certification:

CE-marking? Yes (CE 0120)
Product classification Class IIA (EU)
Quality system certified to ISO 13485
Global Medical Device Nomenclature (GMDN)
 Bacterial/viral filters: 35070
 Heat & moisture exchangers: 46816

Sterilisation:

Sterile? No
Shelf life 5 years

Packaging:

REF*	Shelf box		Shipping carton				MOQ, pcs.
	Qty., pcs.	Dimensions, mm	Qty., pcs.	Dimensions, mm	Gross Wt., kg	Volume, m ³	
761MM	*	*	50	381 x 254 x 203	1.91	0.019	1,250
762MM	*	*	50	381 x 304 x 203	2.27	0.019	1,250
764MM	*	*	50	381 x 254 x 203	2.04	0.019	1,250
864MM	*	*	50	381 x 254 x 203	1.91	0.019	1,250

* = Single packed - no shelf box

* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Storage and handling:



With reservations for misprints and alterations in product specifications. Printed version not controlled.

Manufacturer: **Unomedical Inc.**
 5701-1S, Ware Road
 Address: McAllen, Texas 78503
 USA

European Representative: **Unomedical a/s**
 Address: Kongevejen 2
 3460 Birkerød, Denmark

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.
 The ConvaTec logo is a trademark of ConvaTec Inc.