


All Silicone Foley catheter

Standard, 2-way, hard valve

Intended use:

The all Silicone Foley catheter is a urinary catheter intended for short as well as long term drainage of the bladder. For routine and post-operative drainage. The product is intended for single use only.



List of product range:

All Silicone Foley catheter, Standard, 2-way, hard valve - Regular packaging

CH size	5-10 ml/cc balloon	30 ml/cc balloon
12	MM41151205	
14	MM41151405	
16	MM41151605	MM41151630
18	MM41151805	MM41151830
20	MM41152005	MM41152030
22	MM41152205	MM41152230
24	MM41152405	MM41152430
26	MM41152605	MM41152630

All Silicone Foley catheter, Standard, 2-way, hard valve - Single packed

CH size	1.5 ml /cc pediatric	3 ml/cc pediatric	5-10 ml/cc balloon	30 ml/cc balloon
06	MMS43150603			
08	MMS43150803			
10		MMS43151003		
12			MMS41151205	
14			MMS41151405	
16			MMS41151605	MMS41151630
18			MMS41151805	MMS41151830
20			MMS41152005	MMS41152030
22			MMS41152205	MMS41152230
24			MMS41152405	MMS41152430
26			MMS41152605	MMS41152630

Material specifications:

Contains latex? No
 Contains DEHP? Yes
 Catheter tube and balloon Silicone
 Valve Polypropylene
 Tip Silicone, Barium sulphate
 Sleeve Polyvinyl chloride

Certification:

CE-marking? Yes (CE 0123)
 Product classification Class IIB (EU)
 Quality system certified to ISO 13485
 Global Medical Device Nomenclature (GMDN).34917

Sterilisation:

Sterile? Yes
 Sterilized by Ethylene oxide
 Shelf life 5 years

All Silicone Foley catheter

Standard, 2-way, hard valve

Packaging:

REF*	Shelf box		Shipping carton				MOQ, pcs.
	Qty., pcs	Dimensions, mm	Qty., pcs	Dimensions, mm	Gross Wt., kg	Volume, m ³	
All Silicone Foley catheter, Standard, 2-way, hard valve - Regular packaging							
MM41151205	10	500 x 92 x 52	500	580 x 500 x 527	14.00	0.15	500
MM41151405	10	500 x 92 x 52	500	580 x 500 x 527	14.80	0.15	500
MM41151605	10	500 x 92 x 52	500	580 x 500 x 527	15.60	0.15	500
MM41151630	10	500 x 92 x 52	500	580 x 500 x 527	16.80	0.15	500
MM41151805	10	500 x 92 x 52	500	580 x 500 x 527	17.00	0.15	500
MM41151830	10	500 x 92 x 52	500	580 x 500 x 527	17.00	0.15	500
MM41152005	10	500 x 92 x 52	500	580 x 500 x 527	17.40	0.15	500
MM41152030	10	500 x 92 x 52	500	580 x 500 x 527	18.40	0.15	500
MM41152205	10	500 x 92 x 52	500	580 x 500 x 527	19.00	0.15	500
MM41152230	10	500 x 92 x 52	500	580 x 500 x 527	20.00	0.15	500
MM41152405	10	500 x 92 x 52	500	580 x 500 x 527	22.00	0.15	500
MM41152430	10	500 x 92 x 52	500	580 x 500 x 527	24.00	0.15	500
MM41152605	10	500 x 92 x 52	500	580 x 500 x 527	26.00	0.15	500
MM41152630	10	500 x 92 x 52	500	580 x 500 x 527	17.00	0.15	500

All Silicone Foley catheter, Standard, 2-way, hard valve - Single packed

MMS43150603	*	*	100	576 x 423 x 323	2.50	0.08	500
MMS43150803	*	*	100	577 x 423 x 323	2.50	0.08	500
MMS43151003	*	*	100	578 x 423 x 323	2.50	0.08	500
MMS41151205	*	*	100	576 x 423 x 323	2.50	0.08	500
MMS41151405	*	*	100	577 x 423 x 323	2.50	0.08	500
MMS41151605	*	*	100	578 x 423 x 323	2.50	0.08	500
MMS41151805	*	*	100	579 x 423 x 323	2.50	0.08	500
MMS41152005	*	*	100	580 x 423 x 323	2.50	0.08	500
MMS41152205	*	*	100	581 x 423 x 323	2.50	0.08	500
MMS41152405	*	*	100	582 x 423 x 323	2.50	0.08	500
MMS41152605	*	*	100	583 x 423 x 323	2.50	0.08	500
MMS41151630	*	*	100	576 x 423 x 323	2.50	0.08	500
MMS41151830	*	*	100	577 x 423 x 323	2.50	0.08	500
MMS41152030	*	*	100	578 x 423 x 323	2.50	0.08	500
MMS41152230	*	*	100	579 x 423 x 323	2.50	0.08	500
MMS41152430	*	*	100	580 x 423 x 323	2.50	0.08	500
MMS41152630	*	*	100	581 x 423 x 323	2.50	0.08	500

* = Single packed - no shelf box

* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Storage and handling:



With reservations for misprints and alterations in product specifications. Printed version not controlled.

Manufacturer: **Unomedical Sdn Bhd**
(15290-P)
Address: Bakar Arang Industrial Estate
08000 Sungai Petani, Kedah, Malaysia

European Representative: **Unomedical a/s**
Address: Kongevejen 2
3460 Birkerød, Denmark

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.
The ConvaTec logo is a trademark of ConvaTec Inc.

