


Yankauer Rectoscope

Suction handles and sets

Intended use:

The Yankauer Rectoscope suction handles and sets are for the removal of blood and fluids during rectoscopic surgery. The product is intended for single use only.



List of product range:

	REF*	Connecting tube Shore A 74		Inner/Outer diameter of suction tip mm	Suction handle, length in mm	Vacuum control	Number of relieving eyes	CH/FG (handle)	
		Length in mm	CH/FG						
Yankauer Rectoscope, single bend	3301918x	-	-	4.6/6.5	470	No	0	20	Medium
	3401918x	2,100	24	4.6/6.5	470	No	0	20	Medium
	3402018x	3,500	24	4.6/6.5	470	No	0	20	Medium
	3403618x	2,100	24	4.6/6.5	470	No	0	20	Medium

* Here noted as x, reference codes contain an additional digit (1-5) representing a language code.
Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Material specifications:

Contains latex? No
 Contains DEHP Yes
 Handle Polyvinyl chloride
 Tubing Polyvinyl chloride
 Other components Polyvinyl chloride

Certification:

CE-marking? Yes (CE 0086)
 Product classification Class IIA (EU)
 Quality system certified to ISO 13485
 Global Medical Device
 Nomenclature (GMDN) Handle: 35917/Set: 36778

Sterilisation:

Sterile? Yes
 Sterilized by Ethylene oxide
 Shelf life 5 years

Yankauer Rectoscope

Suction handles and sets

Packaging:

REF*	Shelf box		Shipping carton				MOQ, pcs.
	Qty., pcs	Dimensions, mm	Qty., pcs	Dimensions, mm	Gross Wt., kg	Volume, m ³	
3301918x	50	145 x 124 x 576	300	391 x 304 x 585	6.70	0.072	900
3401918x	*	*	50	596 x 396 x 205	7.80	0.072	150
3402018x	*	*	50	596 x 396 x 205	6.50	0.072	150
3403618x	*	*	50	596 x 396 x 205	7.80	0.072	150

* = Single packed - no shelf box

* Here noted as x, reference codes contain an additional digit (1-5) representing a language code.

Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Storage and handling:



With reservations for misprints and alterations in product specifications. Printed version not controlled.

Manufacturer: **ConvaTec Limited**
First Avenue
Address: Deeside Industrial Park
Deeside, Flintshire, CH5 2NU, UK

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.
The ConvaTec logo is a trademark of ConvaTec Inc.

