



WE CAN

Provide options for more
comfort and flexibility



Introducing

ESTEEM[®] Soft Convex

Flexibility combined with shallow convexity
— designed to deliver comfort



Convex



Flexible,
conforms
to body



Soft




Designed to
minimise
pressure



Designed
to be skin
friendly

ESTEEM® Soft Convex

Esteem®  Soft Convex combines flexibility with a shallow convexity, and is designed to conform well to the body and be comfortable to wear.


Features include:

- **Flexible Convexity:** conforms to the body and helps improve the pouch seal by minimising leaks
- **Easy-View Window:** for accurate positioning and observation
- **Filter:** designed to help minimise odours and ballooning
- **Safe Seal Clipless Closure:** designed to make emptying and cleaning easy*

*Drainable options only.



The Esteem®  Soft Convex system is available in 3 plateau widths to address a range of stoma management issues.



V-Size plateau guidance




Convexity pressure

A narrow plateau width may be used to create pressure directly around the stoma.

Plateau width

A wider plateau width may be used to help flatten out skin creases away from the stoma.

Soft convex base plate depth = 4mm



Esteem® Soft Convex Closed

Pouch size 145 x 205mm

Cut to fit

Product Code	Plateau Size	Cutting Range	Quantity
422360	V1	20-47mm	30
422361	V2	15-40mm	30
422362	V3	10-28mm	30

Pre-cut

Product Code	Plateau Size	Stoma Size	Quantity
422527	V1	25mm	30
422528	V1	30mm	30
422529	V1	35mm	30
422530	V1	40mm	30
422531	V1	45mm	30
422532	V2	20mm	30
422533	V2	25mm	30
422534	V2	30mm	30
422535	V2	35mm	30
422536	V3	15mm	30
422537	V3	20mm	30
422538	V3	25mm	30

Esteem® Soft Convex Drainable

Pouch size 150 x 220mm

Cut to fit

Product Code	Plateau Size	Cutting Range	Quantity
422363	V1	20-47mm	10
422364	V2	15-40mm	10
422365	V3	10-28mm	10

Pre-cut

Product Code	Plateau Size	Stoma Size	Quantity
422539	V1	25mm	10
422540	V1	30mm	10
422541	V1	35mm	10
422542	V1	40mm	10
422543	V1	45mm	10
422544	V2	20mm	10
422545	V2	25mm	10
422546	V2	30mm	10
422547	V2	35mm	10
422548	V3	15mm	10
422549	V3	20mm	10
422550	V3	25mm	10

Esteem® Soft Convex Urostomy

Pouch size 145 x 205mm

Cut to fit

Product Code	Plateau Size	Cutting Range	Quantity
422551	V1	20-47mm	10
422557	V2	15-40mm	10
422562	V3	10-28mm	10

Pre-cut

Product Code	Plateau Size	Stoma Size	Quantity
422552	V1	25mm	10
422553	V1	30mm	10
422554	V1	35mm	10
422555	V1	40mm	10
422556	V1	45mm	10
422558	V2	20mm	10
422559	V2	25mm	10
422560	V2	30mm	10
422561	V2	35mm	10
422563	V3	15mm	10
422564	V3	20mm	10
422565	V3	25mm	10

For Customer Care and to order your FREE samples, call:

Freecall **1800 335 276** (Australia)

0800 441 763 (New Zealand)

or email connection.au@convatec.com

connection.nz@convatec.com

Australia www.convatec.com.au New Zealand www.convatec.co.nz

®/™ indicates trademarks of ConvaTec Inc.

© 2020 ConvaTec Inc. AP-O21239-AU May 2020 O575



Ostomy Care