


Suction Catheters

ProFlo™ and ProFlo™ Plain with Vacutip

Intended use:

Suction catheters are intended for oro-, nasopharyngeal and tracheobronchial suctioning. Must be attached to a jar/canister via tube before being connected to vacuum source. The product is intended for single use only.



List of product range:

	Standard Material (Shore A 78) REF**	Item description	CH	O.D. mm	Length* mm	Color coding
ProFlo™, straight tip with Vacutip, with metric cm-marking on catheter						
	2541005101xx	ProFlo™, straight, Extended Length, Metric	10	3.3	570	Black
	2541006101xx	ProFlo™, straight, Extended Length, Metric	12	4.0	570	White
	2541007101xx	ProFlo™, straight, Extended Length, Metric	14	4.7	570	Green
	2541008101xx	ProFlo™, straight, Extended Length, Metric	16	5.3	570	Orange
ProFlo™ Plain, straight tip with Vacutip, with metric cm-marking on catheter						
	2541001101xx	ProFlo™ Plain, straight, Extended Length, Metric	10	3.3	570	Black
	2541002101xx	ProFlo™ Plain, straight, Extended Length, Metric	12	4.0	570	White
	2541003101xx	ProFlo™ Plain, straight, Extended Length, Metric	14	4.7	570	Green
	2541004101xx	ProFlo™ Plain, straight, Extended Length, Metric	16	5.3	570	Orange

* Effective length

** Here noted as xx, reference codes contain two additional digits (10-50) representing a language code. Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Material specifications:

Contains latex? No
 Contains DEHP? Yes
 Catheter material. Polyvinyl chloride, Styrene-acrylonitrile copolymer

Certification:

CE-marking? Yes (CE 0086)
 Product classification Class IIA (EU)
 Quality system certified to ISO 13485
 Global Medical Device Nomenclature (GMDN)10749

Sterilisation:

Sterile? Yes
 Sterilized by Ethylene oxide
 Shelf life 5 years

Suction Catheters

ProFlo™ and ProFlo™ Plain with Vacutip

Packaging:

Suction Catheters ProFlo™ and ProFlo™ Plain with Vacutip

REF*	Shelf box		Shipping carton				MOQ, pcs.
	Qty., pcs.	Dimensions, mm	Qty., pcs.	Dimensions, mm	Gross Wt., kg	Volume, m ³	
ProFlo™, straight tip with Vacutip							
2541005101xx	100	145 x 124 x 734	600	391 x 304 x 743	9.4	0.088	1,800
2541006101xx	100	145 x 124 x 734	600	391 x 304 x 743	9.4	0.088	1,800
2541007101xx	100	145 x 124 x 734	600	391 x 304 x 743	10.5	0.088	1,800
2541008101xx	100	145 x 124 x 734	600	391 x 304 x 743	10.5	0.088	1,800
ProFlo™ Plain, straight tip with Vacutip							
2541001101xx	100	145 x 124 x 734	600	391 x 304 x 743	8.8	0.088	1,200
2541002101xx	100	145 x 124 x 734	600	391 x 304 x 743	9.4	0.088	1,200
2541003101xx	100	145 x 124 x 734	600	391 x 304 x 743	10.5	0.088	1,200
2541004101xx	100	145 x 124 x 734	600	391 x 304 x 743	11.2	0.088	1,200

* Here noted as xx, reference codes contain two additional digits (10-50) representing a language code.
 Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Storage and handling:



With reservations for misprints and alterations in product specifications. Printed version not controlled.

Manufacturer: **ConvaTec Limited**
 First Avenue
 Address: Deeside Industrial Park
 Deeside, Flintshire, CH5 2NU, UK

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.
 The ConvaTec logo is a trademark of ConvaTec Inc.

